
## LABOROTA 20 eco

Betriebsanleitung  
Instruction Manual  
Mode d'Emploi  
Instrucciones  
de Empleo  
Istruzioni per l'uso

**BEDIENFRONT / FRONT PANEL**  
**PANNEAU FRONTAL DE COMMANDE / FRONTAL DE MANDO**  
**PANNELLO FRONTALE DI COMANDO**

**DEUTSCH**

**Seite 4 - 37**

**ENGLISH**

**page 38 - 71**

**FRANCAISE**

**page 72 - 107**

**ESPAÑOL**

**página 108 - 142**

**ITALIANO**

**página 143 - 179**

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EN 61010-1 + EN 61010-2-010

Thank you for buying a Heidolph Instruments product. This unit has been designed, made and inspected in compliance with DIN EN ISO 61010 for long-term and flawless operation.

## SUMMARY

<b>SUMMARY .....</b>	<b>38</b>
<b>PARTS AND ACCESSORIES SUPPLIED.....</b>	<b>40</b>
<b>GENERAL.....</b>	<b>41</b>
<b>SAFETY INFORMATION .....</b>	<b>42</b>
<b>LABOROTA 20 ECO SPECIFICATIONS.....</b>	<b>43</b>
<b>LABORATORY SET-UP .....</b>	<b>44</b>
<b>SET-UP CHAPTER 1 .....</b>	<b>45</b>
1. Set-up .....	45
2. Install vapor tube and PTFE-seal.....	46
3. Install condensers and expansion vessel .....	47
4. Install collecting flask in cassette.....	48
5. Install inlet pipe with valve and screwed cap .....	49
<b>GLASS-SET - R 15-300-006-21.....</b>	<b>51</b>
<b>GLASS-SET - RC 15-300-006-22 .....</b>	<b>52</b>
<b>GLASS-SET - A 15-300-006-23.....</b>	<b>53</b>
<b>GLASS-SET - AC 15-300-006-24 .....</b>	<b>54</b>
<b>GLASS-SET - A2 15-300-006-25.....</b>	<b>55</b>
<b>GLASS-SET -A2C 15-300-006-26 .....</b>	<b>56</b>
<b>SET-UP CHAPTER 2 .....</b>	<b>57</b>
6. Drive lock.....	57
7. Tightening of Evaporator Flask.....	57
8. Flask supports .....	58
9. Water refill (water bath only).....	58
10. Vacuum and cooling water connection .....	59
11. Fill heating bath .....	60
12. Emptying the heating bath .....	60
13. Overheating protection of the heating bath.....	60
<b>WORKING WITH THE LABOROTA 20 ECO.....</b>	<b>61</b>
<b>OPERATION .....</b>	<b>61</b>
14. Master switch / Stop-Mode switch .....	61
15. Heating bath lift.....	62
16. Speed setting.....	62
17. Heater.....	63
18. Safety functions .....	63
<b>OPTIONS .....</b>	<b>64</b>
19. Vacuum regulator VAC control automatic RS with vacuum valve.....	64
<b>DISTILLATION SETTING .....</b>	<b>65</b>

20. General.....	65
21. Solvent data.....	65
<b>EVAPORATION RATES .....</b>	<b>67</b>
<b>CLEANING &amp; SERVICING.....</b>	<b>68</b>
<b>UNINSTALL, FORWARD &amp; STORE.....</b>	<b>68</b>
<b>DISPOSAL .....</b>	<b>68</b>
<b>TROUBLESHOOTING .....</b>	<b>69</b>
<b>WARRANTY, LIABILITY &amp; COPYRIGHT.....</b>	<b>69</b>
<b>QUESTIONS / REPAIR WORK.....</b>	<b>70</b>
<b>CE-DECLARATION OF CONFORMITY .....</b>	<b>71</b>

**Important information**

**Advise about power cord / mains supply**

**Caution: mandatory action**

**Caution: fire- and explosion hazard**

**Advise about maintenance / repair**

**Caution : Entanglement-rotating parts**

**Caution : Hot surface**

**Caution : keep your hands clear**

## PARTS AND ACCESSORIES SUPPLIED

### Items contained in your carton (s)

Your LABOROTA 20 carton (s) contain the parts and accessories as listed below. Before disposing of the carton, check the contents for all items on the below parts list.

We recommend to keep packing material for future use.

Item	P/N	qty.	for more details refer to page
<b>LABOROTA 20 eco</b>		1	
<b>glass set-R (control)</b>	<b>15-300-006-21</b>	1	<b>page 51</b>
<b>or glass set-RC (control)</b>	<b>15-300-006-22</b>	1	<b>page 52</b>
<b>or glass set-A (control)</b>	<b>15-300-006-23</b>	1	<b>page 52</b>
<b>or glass set-AC (control)</b>	<b>15-300-006-24</b>	1	<b>page 54</b>
<b>or glass set-A2 (control)</b>	<b>15-300-006-25</b>	1	<b>page 55</b>
<b>Large coupling</b>	<b>11-300-003-17</b>	1	<b>page 46, (1)</b>
<b>Support rod Controller compl.</b>	<b>11-300-003-97</b>	1	<b>page 30</b>
<b>Large split ring</b>	<b>23-30-01-02-31</b>	1	<b>page 46, (13)</b>
<b>PTFE O-Ring</b>	<b>23-30-01-01-11</b>	1	<b>page 46, (5)</b>
<b>PTFE-adapter</b>	<b>23-30-01-02-29-1</b>	1	<b>page 46, (4)</b>
<b>PTFE-seal</b>	<b>23-30-01-01-51-0</b>	1	<b>page 46, (6)</b>
<b>Backing ring</b>	<b>23-30-01-02-28</b>	1	<b>page 46, (7)</b>
<b>Small split ring</b>	<b>23-30-01-02-27</b>	1	<b>page 46, (11)</b>
<b>Small coupling</b>	<b>11-300-003-18</b>	1	<b>page 46, (2)</b>
<b>Transportation rod</b>	<b>22-30-01-03-17</b>	2	<b>page 45</b>
<b>Cap</b>	<b>11-300-003-28</b>	4	<b>page 45, (5)</b>
<b>Instruction Manual</b>	<b>01-005-004-68</b>	1	


## Accessories (optional)

Item	P/N
20 l evaporator flask	15-300-003-17
10 l evaporator flask	15-300-003-16
6 l evaporator flask	15-300-003-20
20 l powder flask	15-300-003-19
10 l powder flask	15-300-003-18
Base cart for LABOROTA 20 eco	591-23300-00
Vacuum regulator VAC control automatic RS	591-00341-00
Vacuum valve	591-24000-00
5 liters bath liquid	515-31000-00

### GENERAL



-  Unpack your LR20 carefully.  
Inspect for damage and report such damage or missing parts to your supplier right away.
-  Read your Instruction Manual carefully. Take time to save time while using your LABOROTA 20. Make sure that every user has read and understood the Instruction Manual.
-  Please keep the Instruction Manual in a place easily accessible to every user.

IF ALL ELSE FAILS, READ THESE INSTRUCTIONS !

-  For safety reasons, connect your LABOROTA 20 eco with a grounded outlet only. (also refer to "Lab Installation")

-  **Using de-ionized or distilled water needs to add 0.2 % of Borax ( $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$ ) (corrosion inhibitor)**  
**Please note: The use of tapwater as heating media can lead to calcification which can provoke crevice corrosion. Therefore the user is strongly recommended to clean the heating bath with an adequate polish on a regular basis.**

## SAFETY INFORMATION

-  Please comply with all safety and accident-prevention regulations, as in force for laboratory work !
-  Use extra care when working with rotary evaporators.
-  Use adequate eye protection and protective garments.
-  Use extra care when working with flammable substances; refer to safety data sheets.
-  When connecting your unit with your local power supply, please make sure your unit is designed for your local supply voltage; go by data plate on the unit. Please make sure that the power supply is easily accessible.
-  Turn power switch OFF whenever the unit is not used, or before disconnecting the plug.
-  Repair work is limited to technicians approved by Heidolph Instruments.
-  Use extra care when working in the vicinity of flammable and explosive substances. Motors are of non-sparking type, the unit itself however is not explosion-protected.
-  Caution ! Severe burns may occur when running the heating bath over 50°C. Do not touch hot metal rim of bath, evaporator flask or heating bath fluid.
-  Never run heating bath dry.
-  Your LABOROTA 20 needs a solid level stand.
-  Never attempt to work with damaged or worn-out PTFE gasket (6). Not reaching vacuum pressure wanted may damage your rotary evaporator. Inspect PTFE gasket (6) in regular intervals and replace on condition.
-  Avoid overpressure in the glassware kit while distilling low-boiling solvents. Aerate glassware kit as required. Overpressure in the glassware kit may damage the evaporator or even blow-up parts of the glassware kit.
-  The reliability of the device might be affected, if safety- and installation instructions are being violated.
-  Caution: entanglement-rotating parts
-  Caution : Hot surface
-  Caution: keep your hands clear

## LABOROTA 20 eco SPECIFICATIONS

Supply voltage	3 x 400 V 50 Hz or 3 x 230 V 60 Hz or 1 x 230 V 50/60 Hz
Power rating	4.300 VA
Rotation speed (electronic control)	6 - 160 1/min
Heating bath	400 mm dial., stainless steel V2A
Water bath level regulation	via water feed as per the overflow principle
Temperature range heating bath	25 - 180 ° C
Heating power	4000 watts
Bath lift	motorized
Over temperature & dry running protection	separate safety circuit
Evaporator flask	20 , 10, or 6 liters optional
Collecting flask	10 liters, features bottom drain valve
Pumping capacity required (vacuum pump)	2.5 - 5 m <sup>3</sup> /h (vacuum pump is optional)
Cooling water consumption	150 - 250 l / h (depending on glass set)
Dimensions ( width x depth x height )	970 mm x 620 mm x 1600 mm (with front panel protrusion: 690 mm)
Ambient conditions	0 – 40°C at max. 80% rel. humidity
Weight	approx. 100 kg (does not include glass set)


Illustration including substructure 591-23300-00

## LABORATORY SET-UP

### SAFETY INFORMATION:

**THIS APPLIANCE IS NOT EXPLOSION-PROTECTED. USE EXTREME CARE WHEN USING IN THE VICINITY OF FLAMABLE OR EXPLOSIVE SUBSTANCES.**

**When connecting this product with your local power supply, make sure you got the right thing ! Check VOLTAGE DATA on DATA PLATE (located in the rear).**

This product has to be branched with your local 3-phase network by a skilled electrician.

Wire code:


3 x Phase	GREEN/YELLOW	= protective ground (PE)
	BLUE	= neutral (N)
	BROWN	= phase (L)
	BLACK	= phase (L)
	GREY	= phase (L)
1 x Phase	GREEN	= protective ground (PE)
	WHITE	= neutral (N)
	BLACK	= phase (L)

## SET-UP CHAPTER 1


Carefully unpack all parts, inspect for missing parts by packing list.

### 1. Set-up

- The **LABOROTA 20 eco** with its transportation rods (4) is bolted to the pallet. Remove these screws and carry the unit to its place of installation (CAUTION: 2 persons needed; in this configuration, unit's weight is about 100 kg).
- Place the LABOROTA 20 eco on a stable level stand. In ready-to-go configuration, unit's weight will be about 140 kg ! Ideal stand height above ground is 60 cm, this enables easy access to all controls. We recommend the optional LABOROTA 20 eco base (P/N 591-23300-00 (1)).
- Remove transportation rods (4) and keep for later use.
- Through-holes for carrying rods are plugged with caps (5).
- Install protective cabinet (6) on table (7), opening pointing towards operator. Use screws (8) pre-installed in attaching holes.

## 2. Install vapor tube and PTFE-seal

Unscrew large-size (1) and small (2) coupling rings. Lock driveshaft by turning drive lock (3) CW to ease loosening and tightening of large-size ring (1).

- Remove PTFE-adapter (4) with built-in PTFE-O-ring (5) and PTFE-seal (6)
- Install backing ring (7) and vapor tube (8) in drive head. Put sandwich gasket (9) on vapor tube coupling flange
- Slide small ring (2) over coupling flange of distributor (10). Expand small split ring (11) and slide over distributor (10) coupling flange, shoulder pointing forward
- Secure assembly (10, 2, 11) tight with small coupling ring (2).

Slide PTFE seal (6) on vapor tube (8), sealing lip pointing inward

Never attempt to work with damaged or worn-out PTFE gasket (6). Not reaching vacuum pressure wanted may damage your rotary evaporator. Inspect PTFE gasket (6) in regular intervals and replace on condition.

- Slide PTFE-adapter (4) with built-in PTFE-O-ring (5) into drive module.
- Slide large-size coupling (1) over evaporator flask (12) coupling flange. Expand large split ring (13) and slide over evaporator flask (12) coupling flange, shoulder pointing forward. Set flask support tubes to 10 liter or 20 liter mark and lock with knobs (14); (I = 10 l, II = 20 l)

- Locate evaporator flask on flask lift support tubes (15), tilt flask neck towards drive head, center connecting flange in drive head, so that connecting flange rests against PTFE-O-ring (5) located in PTFE-adaptor (4).


Screw large-size ring (1) tight. When evacuating, double-check for firm seat


- Return flask lift support tubes to initial position (16), CAUTION: basic position uses electronic monitor circuit. Leaving this position will freeze all electric functions of the LABOROTA 20 eco. This situation is reported by an LED on the control panel (flask support).

### 3. Install condensers and expansion vessel

- Loosen thumbscrew (2), lift bracket (1) completely and secure with thumbscrew (2)
- Place rubber side of washer (3) on lower bracket (4), line up holes. Put one more washer (3) with rubber side pointing up on upper neck of condenser or expansion vessel (5)
- Plug lower neck of condenser or expansion vessel (5) through cutout in bracket (4) and washer (3), and tilt to an upright position
- Loosen thumbscrew (2) and slide bracket (1) over upper neck of condenser or expansion vessel. Make vertical adjustment with bracket (4)
- Rotate thumbscrew (2) tight when condenser or expansion vessel are in an exactly upright position.


#### 4. Install collecting flask in cassette

Loosen and remove 4 ea. knurled knobs (1).


- Lift and remove top panel (2).
- Install collecting flask (6) in bottom panel (3). Aerating fitting (4) needs to point forward.
- Reinstall top panel (2) and secure with the 4 knurled knobs (1)
- Install cassette assembly. (with collecting flask) on bracket (8)
- Hold handle (7) of bracket (3), loosen thumbscrew (5), and lower bracket completely; install cassette completely (direction as marked).
- Lift handle (7) until ground joint of collecting flask reaches mating end of condenser, secure with thumbscrew (5)
- use S40 flask clamp to secure both parts in place


## 5. Install inlet pipe with valve and screwed cap

1. Valve sealing point is adjusted by turning an Allen screw clockwise with the Allen wrench (4) supplied.
2. Slide threaded ring (5) and gasket kit (5) over PTFE inlet hose (1), as shown. Secure inlet pipe (2) on distributor (6) with threaded ring (5).
3. Install nipple with (7) on inlet pipe (2) with coupling ring (8). Nipple accepts an ID = 6.6 mm hose.


**Attention: Never install a valve in main return line of cooling water. A closed valve leads to pressure build-up in the glassware which can lead to damage of glass.**

## GLASS-SET - R 15-300-006-21

## components

A	15-300-001-11	condenser, T auto (R) *
B	15-300-002-31	vacuum fitting
C	15-300-002-58	distributor, includes valve (control)
C1	15-300-002-57	valve GL25/165
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	screwed cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve
Z	15-300-001-08	condenser


Use this list when checking your shipment for completeness.

**GLASS-SET - RC 15-300-006-22**

**components**

A	15-300-001-11	condenser, T auto (R) *
B	15-300-002-31	vacuum fitting
C	15-300-002-58	distributor, includes valve (control)
C1	15-300-002-57	valve GL25/165
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
G	15-300-002-59	distributor, includes valve
G1	15-300-002-56	valve GL25/125 lg.
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	screwed cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve


Use this list when checking your shipment for completeness.

**GLASS-SET - A 15-300-006-23**

**components**

A	15-300-001-12	condenser, T auto (A) *
B	15-300-002-31	vacuum fitting
C	15-300-001-09	expansion vessel
D	15-300-002-39	elbow
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
G	15-300-002-40	Tee
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
K	15-300-002-41	distributor, to expansion vessel
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	screwed cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve
Y	15-300-002-60	valve body & valve
Y1	15-300-002-56	valve GL25 / 125 lg.


Use this list when checking your shipment for completeness.

**GLASS-SET - AC 15-300-006-24**

**components**

A	15-300-001-12	condenser, T auto (A) *
B	15-300-002-31	vacuum fitting
C	15-300-001-09	expansion vessel
D	15-300-002-39	elbow
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
G	15-300-002-40	Tee
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
K	15-300-002-41	distributor, to expansion vessel
L	15-300-002-59	distributor, includes valve
L1	15-300-002-56	valve GL25 / 125lg.
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	screwed cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve


Use this list when checking your shipment for completeness.

**GLASS-SET - A2 15-300-006-25**

**components**

A	15-300-001-11	condenser, T auto (R) *
A1	15-300-001-12	condenser, T auto (A2) *
B	15-300-002-31	vacuum fitting
C	15-300-001-09	expansion vessel
D	15-300-002-39	elbow
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
G	15-300-002-40	Tee
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
K	15-300-002-41	distributor, to expansion vessel
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	screwed cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve
Y	15-300-002-60	valve body & valve
Y1	15-300-002-56	valve GL25 / 125 lg.


Use this list when checking your shipment for completeness.

**GLASS-SET -A2C 15-300-006-26**

**components**

A	15-300-001-11	condenser, T auto (R) *
A1	15-300-001-12	condenser, T auto (A2) *
B	15-300-002-31	vacuum fitting
C	15-300-001-09	expansion vessel
D	15-300-002-39	elbow
E	15-300-002-29	vapor tube, ground
F	15-300-003-17	evaporator flask, 20 l
G	15-300-002-40	Tee
H	15-300-004-11	collecting flask, 10 l, features bottom drain
H1	15-300-002-56	valve GL25/125 lg.
I	15-300-002-54	inlet pipe
I1	04-03-01-06-02	hose, PTFE 6.5x0.5
K	15-300-002-41	distributor, to expansion vessel
N	11-300-003-25	GL25 coupling kit
O	15-300-002-45	valve GL18
P	11-300-005-22	nipple for GL14
Q	23-09-03-01-27	coupling ring GL14
R	23-30-01-02-64	nipple for GL18
S	23-09-03-01-23	coupling ring GL18
T	23-09-03-01-24	threaded cap GL18
U	11-300-002-94	threaded joint KS 64, coated
V	11-300-002-95	flask clamp for S40/S41
W	23-30-01-01-09	sandwich gasket 50 / 69 / 3
X	23-30-01-03-31	spacer sleeve
Y	15-300-002-59	distributor, includes valve
Y1	15-300-002-56	alve GL25 / 125 lg.


Use this list when checking your shipment for completeness.

## SET-UP CHAPTER 2

## 6. Drive lock


The drive lock facilitates opening and closing the large coupling ring (3), holding the flask (4).

- To engage lock: turn knob (5) CW until **“lock”** appears. Manually rotate drive head, until lock snaps into place. **Rotation Stop LED (6)** on front panel illuminates.
- To release lock: turn knob (5) CCW until **“unlock”** appears. **Rotation Stop LED (6)** on front panel dims.


## 7. Tightening of Evaporator Flask

Evaporator Flask must be sufficiently tightened with large coupling ring (3).

**Check that the screws are fixed firmly; if the cap nuts are not fastened tightly enough, this will cause leakage and functional disturbances!**

## 8. Flask supports

The flask supports make flask (1) removal (even filled) an easy job for one person.

3. **Uninstall flask:** Engage drive lock. Slide support pipes (4) to flask (1) from both sides. Check for firm contact and lock support pipes in place by tightening thumbscrews (5) CW. Unscrew large coupling ring (3, not shown) by turning CCW. Tilt flask neck with large coupling ring (3) forward. Remove flask carefully (coupling ring to remain on flask).

4. **Install flask:** Engage drive lock. Move support pipes (4) (both sides) to the respective flask size mark (center of thumbscrew = mark). Mark I = 10 l flask;

mark II = 20 l flask. Lock support pipes in place by tightening thumbscrews (5) CW. Put flask (1), holding large coupling ring (3), on support pipes carefully. Swing flask axis with coupling ring in-line with axis of rotation, until flange is centered in drive head. Secure PTFE-adaptor (not shown) with PTFE-O-ring tight by CW turn. Double-check for firm seat. Return support pipes to initial position (7). Release drive lock. **Caution:** Basic position of support pipes is monitored by the electronic circuit; when leaving this position, all electric functions of the system will freeze.

This situation is reported by the **Flask Support** LED on the control panel.

**Before starting the rotation, disengage the blocking device.**

## 22. Water refill (water bath only)


The LABOROTA 20 eco features a water refill system.

The water feed consists of 2 components:

- Overflow tube with drain device (10)
- Water feed valve with water connection and feed tube (11)

Fasten a suitable hose (dia. = 8 mm, temperature-resistant up to 100°C) securely to the hose shaft of the drain device (12) using a hose clamp. Place the end of the hose in a sufficiently large catch tank (13) or water overflow basin. In so doing, make absolutely sure that the end of the hose is at a lower level than the hose shaft so that excess water can flow off.

Connect the water feed valve to the water supply line using a pressure resistant water hose, and secure the hose at the hose shaft (14) using hose clamp.


At the beginning of distillation open the water feed valve (11) just enough to replace the evaporated water. Any excess water is conducted off through the overflow tube (10) into the catch tank (13).

 **Caution:** To avoid wasting water, open the water feed valve during distillation and close it again after distillation.

**23. Vacuum and cooling water connection**

 **The cooling coils of the cooler may only be exposed to a maximum cooling water pressure of 2 bar. So that higher pressure can never develop, no valve and no restrictor may be installed in the cooling water return.**

The necessary vacuum and cooling water hoses are attached to the cooler using the hose couplings (15) and screw caps (16). It is mandatory to secure the cooling water hoses (17) with hose clamps (18)!

## 24. Fill heating bath

The heating bath can either be filled with water or a suitable heat transfer oil.

Using de-ionized or distilled water needs to add 0.2 % of Borax ( $\text{Na}_2\text{B}_4\text{O}_7 \cdot 10\text{H}_2\text{O}$ ) (corrosion inhibitor)

**Please note:** The use of tapwater as heating media can lead to calcification which can provoke crevice corrosion. Therefore the user is strongly recommended to clean the heating bath with an adequate polish on a regular basis.

The filling amount must be checked during the filling procedure. When the flask is mounted and the heating bath lift is in the upper position, the ideal filling level is approx. 3 cm below the edge of the heating bath. Attention must always be given to the volume displacement of the evaporating flask.


## 25. Emptying the heating bath

**Caution:** The bath fluid must be cooled down to room temperature before emptying, otherwise there is a danger of the hot bath fluid scalding and damaging the material.

On the LABOROTA 20 eco there is a drain cock (1) on the back of the device. The valve is closed when the black lever is at 90° to the direction of flow.

This valve drains bath fill into a tank (not supplied).


## 26. Overheating protection of the heating bath

If the overheating protection of the heating bath is actuated (LED **Overheat Prot.** illuminates), the bath must be cooled down and the cause eliminated.

Possible cause: Lack of heating bath fluid.

Descend heating bath lift to LOW.

Turn master switch OFF.

Reset over temperature protector by introducing an insulated screwdriver through the little hole at the lower right end of the unit. (remove black plastic plug to gain access through panel)

## WORKING WITH THE LABOROTA 20 ECO

Use rotary evaporators with appropriate care. Comply with standard laboratory safety practices.

Be cautious when working with flammable solvents; motors and limit switches are of non-sparking type but we cannot guarantee that the unit is explosion-proof.

### SAFETY INFORMATION:


**WHEN CONNECTING YOUR UNIT WITH YOUR LOCAL POWER SUPPLY, PLEASE MAKE SURE YOUR UNIT IS DESIGNED FOR YOUR LOCAL SUPPLY VOLTAGE; GO BY DATA PLATE ON THE UNIT.**

**CONNECT YOUR ROTARY EVAPORATOR WITH A GROUNDED OUTLET ONLY.**

## OPERATION

### 27. Master switch / Stop-Mode switch

The master switch disconnects all wires from power supply, and serves as a RESET button. (after **Stop-Mode**).

The **Stop-Mode** switch discontinues operation of the rotary evaporator immediately.

Stop rotation

Lowers heating bath

Opens aerating valve (vacuum valve closes)

Turns OFF heating circuit

**Stop-Mode** LED turns ON

System re-start after an **Stop-Mode** situation requires system RESET. Make sure the problem is fixed before resetting. See an expert, if required.


To Reset, turn **Master switch** OFF & ON.

## 28. Heating bath lift

Heating bath lift RAISES, as long as the **Lift UP** key is held depressed.

Heating bath lift LOWERS, as long as the **Lift DOWN** key is held depressed.

At upper and lower limits, lift motion stops automatically.

## 29. Speed setting

The rotation is set using the "Speed" knob. On the display "Speed rpm" the speed is shown in rpm (revolutions per minute).

By pressing the **Start Stop button** the rotation is started. The LED next to the Start Stop button indicates that it is switched on. Press once again to stop the rotation.

**CAUTION:** for safety reasons, start rotation at low speeds. Rotational speed can be increased during operation.

Never exceed a speed of 100 rpm, while working with a 20 l evaporation flask.

If rotation does not start with the START/STOP key, please check for illumination of one or more of the Safety Cover, Flask Support or Rotation Stop LED's.

If this occurs check if the enclosure is open, flask support not in one of its extreme positions, or the drive drive lock engaged.

For more information refer to Troubleshooting Section.

## 17. Heater

The heating bath temperature is set using the "Bath Temp" knob. On the display "Bath Temp °C" the temperature is given in °C. It always indicates the actual temperature of the heating bath fluid. When setting the required temperature, the display jumps to the required temperature indicator which then begins to flash.


**Caution:** The heating bath has a temperature range of 20 – 180°C. Make sure that there is oil in the heating bath container when setting to temperatures of over 100°C. As a warning sign, the information "OIL" appears on the display for 2 seconds when 100°C are exceeded

The heating is switched on by pressing the **Start Stop button**. The LED next to the Start Stop button indicates that it is switched on. Press once again to stop the heating.

## 18. Safety functions

All of the following cause LED ON situations.

- **Flask Support LED**: Flask support is not in its extreme position (either side). Rotation won't start, or rotation is going to stop, lift won't work. Move flask support to extreme position (both sides) and secure with knurled knobs, (page 58 (5)).
- **Rotation Stop LED**: Drive lock engaged, rotation can't start. Release lock by turning knob CCW.
- If the **LED Overheat Prot.** illuminates, the overheating protection for the heating is actuated. The Stop - Mode function is initiated. (Heating bath lift travels downwards, heating and rotation come to a halt).  
Heater won't start, heater is disconnected.  
Possible reasons: heating bath level insufficient, heating bath dry; refill. Heater is broken. After eliminating the fault, allow the bath to cool down and push in the overheating protection with a screwdriver. (see also page 26)
- If the **LED Stop - Mode** illuminates, the Stop - Mode switch was pressed or the overheating protection was actuated (see previous paragraph).  
Heating bath lift descends automatically;  
Heater and rotation OFF;  
Reset with master switch OFF and ON



## OPTIONS

**19. Vacuum regulator VAC control automatic RS with vacuum valve**

For regulating the vacuum, the vacuum regulator no. 591 00341 00 and the vacuum valve 591 24000 00 can be supplied.

The following figure shows how the vacuum regulator VAC control automatic (20) and the vacuum valve (21) can be attached to the Laborota 20 eco.

## DISTILLATION SETTING

### 20. General

Using rotational evaporators for thermal breakdown of substances needs some basic considerations about optimizing the distillation process. Parameters stipulated below are to be considered for making optimal settings.

- **Rotational speed, evaporator flask**

Increasing rotational speed will increase rate of distillation. Shorter distillation time preserves thermal stability of your substances to be distilled.

- **Difference between heating bath and vapor temperatures:**

This temperature difference should not be less than 20° C to maintain proper rate of distillation. As long as your substances in distillation allow, you may even chose a higher temperature difference.

Rule of thumb: doubling actual temperature difference doubles rate of distillation.

- **Temperature difference between vapor temperature and coolant:**

To maintain adequate condensation rates, we recommend a 20° C temperature difference between the solvent vapor temperature and coolant. A low temperature difference will affect solvent recovery. For substances with high evaporation heat, you should increase temperature difference.

### 21. Solvent data

- Examples for using table and chart:

- The table displays essential substance data, the chart makes reference to the relationship between pressure and boiling point of some solvents. In this case, 1/T is plotted as a function of log p in accordance with the Clausius-Clapeyron equation.

**a) Distillation of a solvent contained in the table:**

- On the vacuum controller that can be supplied as an accessory, the value given in the column "Vacuum for Spd at 40°C" is set as the required value for Set Vacuum.

**b) Distillation of a solvent contained in the chart, with a boiling point other than 40°C:**

- Mark the boiling point desired on the temperature axis of the chart
- Draw a line to the right to get an intersection with the solvent
- Vacuum pressure required for this one situation is found by drawing a vertical line from this intersection down to the horizontal axis.

**c) Distillation of a solvent not contained in the table or chart:**


- The following notes may be helpful to determine the appropriate vacuum pressure:
- The line gradient is a function of enthalpy of evaporation. It is similar for chemically related substances with boiling points close to each other. Hence, the plotted lines may be used as a guide for substances with slightly different boiling points.
- Using a water-jet or diaphragm-type pump, you may lower boiling point by approx. 100°C.
- Rule of thumb: decrease pressure to a half will lower boiling point by approx. 15°C.

Avoid overpressure in the glassware kit while distilling low-boiling solvents. Aerate glassware kit as required. Overpressure in the glassware kit may damage the evaporator or even blow-up parts of the glassware kit.


Solvent	Empirical formula	MW [g/mol]	Boil. [°C]	$\Delta H_{\text{vap}}$ [J/g]	Vacuum for boil.at 40 °C	
					[mbar]	[mm(Hg)]
Acetone	C <sub>3</sub> H <sub>6</sub> O	58.08	56.5	550	556	387
Acetonitrile	C <sub>2</sub> H <sub>3</sub> N	41.05	81.8	833	230	173
Benzene	C <sub>6</sub> H <sub>6</sub>	78.11	80.1	549	236	177
n-butanol (butyl alcohol)	C <sub>4</sub> H <sub>10</sub> O	74.12	117.5	619	25	19
Tert. butanol (tert. butyl alcohol)	C <sub>4</sub> H <sub>10</sub> O	74.12	82.9	588	130	98
2-butanone (methyl ethyl ketone)	C <sub>4</sub> H <sub>8</sub> O	72.11	79.6	473	243	182
Chlorobenzene	C <sub>6</sub> H <sub>5</sub> Cl	112.60	132.2	375	36	27
Cyclohexane	C <sub>6</sub> H <sub>12</sub>	84.16	80.7	389	235	176
1,2-dichloroethane	C <sub>2</sub> H <sub>4</sub> Cl <sub>2</sub>	98.96	82.4	336	210	158
1,2-dichloroethylene (cis)	C <sub>2</sub> H <sub>2</sub> Cl <sub>2</sub>	96.94	59.0	320	479	134
1,2-dichloroethylene (trans)	C <sub>2</sub> H <sub>2</sub> Cl <sub>2</sub>	96.94	47.8	313	751	563
Dichloromethane (methylene dichloride)	CH <sub>2</sub> Cl <sub>2</sub>	84.93	40.7	373	atm.	atm.
Diethyl ether	C <sub>4</sub> H <sub>10</sub> O	74.12	34.6	392	atm.	atm.
Diisopropyl ether	C <sub>6</sub> H <sub>14</sub> O	102.20	67.5	318	375	281
Dimethyl formamide	C <sub>3</sub> H <sub>7</sub> NO	73.09	153.0		11	8
1,4-dioxan	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88.11	101.1	406	107	80
Ethanol	C <sub>2</sub> H <sub>6</sub> O	46.07	78.4	879	175	131
Ethyl acetate	C <sub>4</sub> H <sub>8</sub> O <sub>2</sub>	88.11	77.1	394	240	180
Heptane	C <sub>7</sub> H <sub>16</sub>	85.09	98.4	439	120	90
Hexane	C <sub>6</sub> H <sub>14</sub>	86.18	68.7	370	335	251
Methanol	CH <sub>4</sub> O	32.04	64.7	1225	337	253
3-methyl-1-butanol (isoamyl alcohol)	C <sub>5</sub> H <sub>12</sub> O	88.15	130.6	593	14	11
Pentachloroethane	C <sub>2</sub> HCl <sub>5</sub>	202.30	160.5	203	13	10
Pentane	C <sub>5</sub> H <sub>12</sub>	72.15	36.1	382	atm.	atm.
n-pentanol (amyl alcohol)	C <sub>5</sub> H <sub>12</sub> O	88.15	137.8	593	11	8
1-propanol (n-propyl alcohol)	C <sub>3</sub> H <sub>8</sub> O	60.10	97.8	787	67	50
2-propanol (isopropyl alcohol)	C <sub>3</sub> H <sub>8</sub> O	60.10	82.5	701	137	103
1,1,2,2-tetrachloroethane	C <sub>2</sub> H <sub>2</sub> Cl <sub>4</sub>	167.90	145.9	247	35	26
Tetrachloroethylene	C <sub>2</sub> Cl <sub>4</sub>	165.80	120.8	233	53	40
Tetrachloromethane (carbon tetrachloride)	CCl <sub>4</sub>	153.80	76.7	225	271	203
Tetrahydrofuran	C <sub>4</sub> H <sub>8</sub> O	72.11	66.0	–	357	268
Toluene	C <sub>7</sub> H <sub>8</sub>	92.14	110.6	425	77	58
1,1,1-trichloroethane	C <sub>2</sub> H <sub>3</sub> Cl <sub>3</sub>	133.40	74.1	251	300	225
Trichloroethylene	C <sub>2</sub> HCl <sub>3</sub>	131.40	86.7	265	183	137
Trichloromethane (chloroform)	CHCl <sub>3</sub>	119.40	61.3	263	474	356
Water	H <sub>2</sub> O	18.02	100.0	2259	72	54
Xylene (isomer mixture)	C <sub>8</sub> H <sub>10</sub>	106.20	137–143	390	25	19

Conversion factor [mbar] to [mm(Hg)]:  $[mm(Hg)] \approx \frac{3}{4}[mbar]$

**EVAPORATION RATES**

Evaporation rates Laborota 20

## CLEANING & SERVICING

**Cleaning:** wipe housing clean with a damp cloth (add some mild liquid soap).

### Note

To avoid damage to the surface finish, avoid using chlorine bleach, chlorine-based detergents, abrasive substances, ammonia, rags or cleaning agents containing metal particles.

**Rust stains appearing on the bottom of the hot bath after long-term use are caused by foreign matter deposits (iron particles or rust particles). Remove such rust stains and calcareous deposits with a conventional phosphoric acid bathroom cleaner (decalcifying agent).**

**Never attempt to work with damaged or worn-out PTFE gasket (refer to page 46 (6)). Not reaching vacuum pressure wanted may damage your rotary evaporator. Inspect PTFE gasket (6) in regular intervals and replace on condition.**

The unit is maintenance-free. Repair work is limited to technicians approved or appointed by Heidolph Instruments. Please call your local Heidolph Instruments Dealer or a Heidolph Instruments Field Representative (also refer to page 70).

## UNINSTALL, FORWARD & STORE

### Uninstall

Please switch off the device and disconnect it from the power supply

### Forward & Store

1. We recommend to store the item and its components in its original box, or a similar container that offers adequate protection against damage in transit. Tape the box securely.
2. Store the item in a dry place.

### Caution

Do not jolt or shake the item during transport.

## DISPOSAL

For disposal, please comply with your local or national regulations.  
Separate by metal, plastic, etc.  
Packing material to be treated as described above (material split).

## TROUBLESHOOTING

- Master switch won't light when turning ON  
check power cord  
master switch broken (call service)
- Poor vacuum pressure  
inspect gaskets and fittings  
check ground fittings, coat with vacuum grease, as required  
PTFE seal installed the wrong way. Sealing lip needs to point towards drive head  
check line to vacuum pump
- Heating bath lift won't travel up  
also refer to safety functions section  
lift mechanism broken  
upper limit switch broken
- Heating bath lift won't travel down  
lift mechanism broken  
lower limit switch broken
- No rotation  
also refer to safety functions section  
motor broken
- Fault message "Er1" on the display Bath Temp: Sensor short circuit, temperature sensor faulty
- Fault message "Er2" on the display Bath Temp: Sensor broken, temperature sensor faulty
- **LED Stop-Mode, LED Flask Support, LED Rotation Stop or LED Overheat Prot.**  
Illuminates, refer to safety functions (page 63)

Whenever facing a problem that cannot be fixed by the above information, call your authorized Heidolph Instruments dealer right away.

## WARRANTY, LIABILITY & COPYRIGHT

### Warranty

Heidolph Instruments provides a three-year warranty on the products described here (with the exception of glass and consumable parts) if registered with enclosed warranty card or via internet ([www.heidolph.com](http://www.heidolph.com)). Warranty starts with the date of registration. Without registration warranty starts according to serial number. This warranty covers defects in materials and workmanship.

Transit damage is excluded from this warranty.

To obtain such warranty service, contact Heidolph Instruments (phone: +49 – 9122 - 9920-68) or your local Heidolph Instruments Dealer. If defects in material or workmanship are found, your item will be repaired or replaced at no charge.

Misuse, abuse, neglect or improper installation are not covered by this warranty promise. Alterations to the present warranty promise need Heidolph Instruments' consent in writing.

### Exclusion of liability

Heidolph Instruments cannot be held liable for damage from improper use or misuse. Remedy for consequential damage is excluded.

### Copyright

Heidolph Instruments is the copyright holder for all texts and images in this manual.

## QUESTIONS / REPAIR WORK

If any **aspect** of installation, operation or maintenance remains unanswered in the present Manual, please get in touch with the following address:

For repairs please call Heidolph Instruments (phone: ++49-9122-9920-68) or your local, authorized Heidolph Instruments Dealer.

### Note

You will receive approval for sending your defective item to the following address:

**Heidolph Instruments GmbH & Co. KG**  
**Lab Equipment Sales**  
**Walpersdorfer Str. 12**  
**D-91126 Schwabach / Germany**  
**phone: ++49-9122-9920-68 Fax: ++49-9122-9920-65**  
**E-Mail: Heidolph.Instruments@Heidolph.de**

### Safety Information

When shipping items for repair that may have been contaminated by hazardous substances, please:

- advise exact *substance*
- take adequate *protective* action towards our parts receiving and service personnel
- *mark* the pack IAW Hazardous Materials Act

## CE-DECLARATION OF CONFORMITY

We herewith declare that the present product complies with the following standards and harmonized documents:

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**© Heidolph Instruments GmbH & Co. KG**