# Intelligent Evaporation

Leading Safety Standards

Superior Ease of Use

Reduced Cost of Ownership

## Operating Manual

Page 2

Rotary Evaporator Hei-VAP Value, Hei-VAP Value Digital, Hei-VAP Advantage

Operating manual must be read before initial start-up.

Please follow the safety instructions provided.

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<b>About this Document</b>	<b>4</b>
Version and variants .....	4
About this manual .....	4
Symbols and keywords used .....	5
<b>Safety Instructions</b>	<b>6</b>
General safety instructions / Consignes de sécurité générales.....	6
EU Declaration of Conformity / Déclaration de conformité UE .....	6
Intended use / Utilisation conforme.....	7
Installation/electrical safety / Installation, sécurité électrique .....	7
Personnel qualification / Qualification du personnel.....	8
Operating company's obligations / Obligations de la société exploitante .....	8
Installation site / Lieu de montage .....	8
Modifications to the device / Modifications de l'appareil.....	9
Personnel safety / Sécurité du personnel.....	9
Safety during use / Sécurité durant l'utilisation.....	10
Disposal / Mise au rebut .....	11
<b>Device Description</b>	<b>12</b>
Device overview .....	12
Hei-VAP Value, Hei-VAP Advantage: Overall view.....	12
Hei-VAP Value: Control panel.....	12
Hei-VAP Value Digital: Control panel.....	13
Hei-VAP Advantage: Control panel .....	13
<b>Start-Up</b>	<b>14</b>
Connecting/disconnecting power cord.....	14
Switching the device on and off.....	14
Heating bath and lift .....	15
Fill heating bath .....	16
Connect heating bath .....	16
Empty heating bath .....	17
Menu navigation.....	18
<b>Operation</b>	<b>19</b>
Adjusting flask rotation .....	19
Setting bath temperature .....	21
Determine vapor temperature .....	23
Setting duration for process .....	23
Evaporating and receiving flask .....	24
Mount and remove evaporating flask .....	25
Mount and remove receiving flask.....	26
Dispensing sample.....	26
Ventilation .....	26

<b>Cleaning and Maintenance</b>	<b>27</b>
Cleaning.....	27
Maintenance .....	28
<b>Troubleshooting</b>	<b>29</b>
Malfunctions and troubleshooting.....	29
Overheat protection: reset.....	31
<b>Assembly</b>	<b>32</b>
Electrical connections.....	32
Connexions Electriques.....	33
Mounting device .....	34
Assembly parts.....	34
Transportation lock .....	35
Heating bath .....	36
Vapor tube.....	36
Condenser .....	38
Inlet tubes.....	43
Vapor temperature sensor, AUTOaccurate sensor.....	44
Evaporating flask angle and depth .....	44
Control panel .....	46
Connecting peripheral devices .....	47
Connect vacuum.....	47
Connect cooling water .....	49
Mount tubes.....	50
Additional accessories .....	50
<b>Disassembly and Storage</b>	<b>51</b>
Dismantling, transportation and storage .....	51
Dismantling .....	51
Transportation and storage .....	52
<b>Accessories and Spare Parts</b>	<b>53</b>
Scope of delivery.....	53
Accessories.....	53
<b>Attachments</b>	<b>54</b>
Technical data .....	54
<b>Service</b>	<b>55</b>
Contact/Technical Support .....	55
Warranty .....	56
Confirmation of condition.....	57
Déclaration de non-opposition .....	58

## ➤ Version and variants

### Version

This manual describes the function, operation and maintenance of the Hei-VAP Value, Hei-VAP Value Digital and Hei-VAP Advantage evaporators.

Version	Alteration date
1.2	01/2018

### Variants

The devices are available in different configurations. Certain features and functions apply to specific configurations only. The varieties are documented in this manual.

## ➤ About this manual

This operating manual is a component of the device described here.

- Please read this manual carefully and obey all safety and warning notices.
- Ensure that every operator reads this manual.
- Ensure that this manual is accessible for every operator.
- Pass on the operating manual to the subsequent owner.

For the current version of this operating manual in pdf format, please go to:

<http://www.heidolph.com/Support/Operation Manuals/Rotary Evaporators/Hei-VAP Series/Hei-VAP Value, Hei-VAP Value Digital and Hei-VAP Advantage.>

In addition, please observe the regional regulations.

### ➤ Symbols and keywords used

In this manual you will find standardized terms and signs to warn you of possible dangers or give important hints. Please take special notice to these signs to avoid personal injury or damage to the appliance.

The following terms and type of signs (basic drawing elements) are used:

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This device complies with the following EC-Directives\*:

- 2006/42/EC Machinery Directive
- 2014/30/EU Electromagnetic Compatibility Directive

\*Also see attachment „EU Declaration of Conformity“.

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The device has been constructed according to state-of-the-art technology and recognized safety regulations. However, risks may still arise during installation, operation and maintenance./

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### **The device may only be used under the following circumstances:/**

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## **Personnel qualification / Qualification du personnel**

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## Safety Instructions

- The device may only be operated in enclosed spaces and under the following environmental conditions: /

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## Modifications to the device / Modifications de l'appareil

- No unauthorized changes may be made to the device.
- No parts may be used which have not been approved by the manufacturer.
- Unauthorized changes result in the EC Declaration of Conformity losing its validity, and the appliance may no longer be operated.
- The manufacturer is not liable for any damage, danger or injuries that result from unauthorized changes or from operating the device other than as described in this manual.

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## ➤ Device overview

### Hei-VAP Value, Hei-VAP Advantage: Overall view

### Hei-VAP Value: Control panel


### Hei-VAP Value Digital: Control panel

### Hei-VAP Advantage: Control panel

### ➤ Connecting/disconnecting power cord

#### Connect power cord

The power cord comes with a three-wire plug and a recess on the bottom side.

Unit plug

Power cord

The socket is located on the back side of the device.

- ✓ The device is switched OFF.
- Connect plug to the socket with the recess facing down.
- Push plug up against socket until securely attached.
- Connect plug to power socket.

Unit socket

#### Disconnect power cord

- Disconnect plug from power socket first.
- Then disconnect plug from socket on device.

### ➤ Switching the device on and off

The ON/OFF switch is situated on the left side of the device beneath the control panel.

#### Switch device on


- Press the ON/OFF switch on the left side.

##### Hei-VAP Value

- ✓ The operating indicator light switches on.

##### Hei-VAP Value Digital, Hei-VAP Advantage

- ✓ The display switches on and shows set parameters.
- ✓ The device is switched on.

ON/OFF switch: ON

#### Switch device off

- Press the ON/OFF switch on the left side.

##### Hei-VAP Value

- ✓ The operating indicator light switches off.

##### Hei-VAP Value Digital, Hei-VAP Advantage

- ✓ The display switches off.
- ✓ The device powers off.

### ➤ Heating bath and lift

You can fill the heating bath with various heating bath fluids:

- Tap water (no deionized or distilled water)
- Water soluble polyethylene glycol
- Oil with low viscosity (40 cP) and flashpoint > 285°C.

#### **Danger/ Danger Risk of explosion! / Danger d'explosion !**

Heating bath fluid oil with a flashpoint < 285 °C might ignite and lead to uncontrollable thermic reactions and explosion./

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The marks inside heating the bath show the minimum and maximum level for heating bath fluid with evaporating flask completely immersed.

Evaporating flask is immersed by moving up and down via the heating bath lift.

When heating bath lift is moved downwards the evaporating flask displaces already filled in fluid.

To prevent heating bath from overflow:

- Move down the lift with evaporating flask mounted before filling bath.
- Fill in heating bath fluid till the maximum level mark is reached.

Do not let the evaporating flask touch the bottom or the edge of the heating bath.

(To adjust the evaporating flask depth and flask angle please view chapter "Assembly", section "Evaporating flask: Adjust Evaporating flask angle and depth".)

## Start-Up

### Immerse Flask: Hand lift


- Push down hand lift lever.
- Move pushed hand lift lever to the left to lower or to the right to raise lift.
- Release hand lift at the desired position.
- ✓ The hand lift locks into place.

Hand lift lever: upper lift position

### Immerse Flask: Motor lift

- ✓ Device is switched on.
- Push lift key down  or up  to move motor lift.
- Release lift key at desired position.
- ✓ The motor lift stops.

Lift keys: up, down

## Fill heating bath

- ✓ Evaporating flask is immersed.
- Fill in heating bath fluid till maximum level mark is reached.
- ✓ The level of fluid is between the minimum and the maximum mark.

## Connect heating bath

### Warning:

#### Risk of crushing! Risk of scalding!

When device is switched on, rotation or heating can be started accidentally and result in severe injuries.

Only connect heating bath to base device when the device is switched off.

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## Start-Up

The heating bath plug is located on the right side in the back of the heating bath.


- ✓ Device is switched off.
- ✓ Heating bath is locked into place at the base of the device.

- Connect 7-pole socket to plug of the heating bath with the recess facing up.
- Fix with knurled nut turning it clockwise.

Knurled nut

Plug heating bath

Recess

Feather

## Empty heating bath

### Warning:

#### Risk of scalding! Risk of slipping!

When emptying the heating bath, you may burn yourself or slip on spilled fluid.

Before emptying heating bath:

- Remove evaporating flask.
- Disconnect base device and heating bath.
- Wait until the heating bath fluid cools down to <math><50^{\circ}\text{C}</math>.
- Please wear protective clothes.
- Lift the heating bath only by the insulated handle./

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- ✓ Evaporating flask is removed.
- Turn knurled nut counter clockwise to unlock it.
- Remove socket from heating bath.
- Empty heating bath.
- Lift bath carefully at both handles and empty.


## ➤ Menu navigation

### Control panel Hei-VAP Value Digital

Start screen shows the last set and the actual temperature.

- ✓ Both rotary knobs are turned left as far as they will go.
- ✓ The device is switched on.
- ✓ The display shows set and actual temperature.
- ➔ Turn rotation knob or temperature knob clockwise.
- ✓ When adjusting temperature the display shows the setting.
- ➔ Press ON/OFF key above knob.
- ✓ ON/OFF key lights up and set process starts.
- ➔ Press ON/OFF key again.
- ✓ Light of ON/OFF key switches off.
- ✓ The process is discontinued.

### Control panel Hei-VAP Advantage

The start screen shows the actual values. Adjust parameters temperature and rotation with select knob Hei-GUIDE.

- ➔ Turn select knob Hei-GUIDE clockwise slowly to the right.
- ✓ Selected parameter is framed.
- ✓ The frame jumps to next parameter when turning select knob.
- ➔ Press select knob Hei-GUIDE.
- ✓ The addressed parameter field is shown.
- ➔ Turn select knob Hei-GUIDE to adjust parameter.
- ➔ Press select knob to confirm setting.
- ➔ Turn select knob Hei-GUIDE to select and adjust next parameter.
- ➔ Start heating and rotation by pushing appropriate ON/OFF key.
- ✓ When heating is switched on the heating key lights up.
- ✓ A preset timer starts as soon as rotation is started.

## ➤ Adjusting flask rotation

By increasing flask rotation you can increase your evaporation rate. A shorter evaporation time may prevent thermal damage to your sample.

Depending on flask size and heating bath temperature different rotation speeds may be suitable.

### Speed range and setting

Rotation speed can be set from 10 rpm - 280 rpm.

#### **Warning:** **Risk of scalding and slipping!**

At high flask rotation speeds bath fluid may splash or spill over and result in a slippery floor and/or severe scalding.

If you are operating the device without a safety hood:

- Run a test with your heating bath fluid temperature and rotation speed and document any size of splashing.
- Ensure that this documentation is communicated and accessible to all individuals.
- Wear safety clothes if you are working in near proximity to the device.
- Pay attention to any amount of heating bath fluid on the floor.

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→ Test your selected fluid

- with slowly increasing temperatures
- with different evaporating flasks
- with increasing rotation speeds

→ Note these results and make them accessible for all contributors.

Depending on the oil and age of oil you are utilizing the heating bath temperature may change the viscosity and the volume of the bath fluid.

Keep records of the exact description/manufacturer information and period of use and repeat the documentation in the determined time interval.

# Operation

## Hei-VAP Value

- ✓ Rotation knob is turned to the left as far as it will go.
- ✓ Device is switched on.
- ✓ Evaporating flask is immersed into the heating bath.
- Press rotation key .
- Enter desired rpm with rotation knob.
- Press rotation key to stop the rotation.

Rotation key

Rotation knob

## Hei-VAP Advantage

- ✓ Device is switched on.
- Select „rotation“ field with select knob Hei-GUIDE.
- Enter desired rpm and confirm.
- Press rotation key  to start rotation.
- Adjust speed via field „rotation“.
- ✓ Actual rpm is shown in the field „rotation“.
- ✓ Time of your process is shown under „timer“ in minutes. An up arrow is displayed next to the time.
- Press rotation key  to stop rotation.

„rotation“ field

Select knob Hei-GUIDE

Rotation key

## ➤ Setting bath temperature

### Possible temperature range:

- Water bath application 20°C - 100°C
- Oil bath application 20°C - 210°C

With a pre-set temperature above 75 °C to 80 °C we recommend oil or polyethylen glycol as heating bath fluid.

Set "Oil" as medium if you work with polyethylen glycol. Please follow temperature specification provided by manufacturer.

### **Danger/ Danger Risk of explosion! / Danger d'explosion !**

Potential loss of life/ Danger de mort

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### **Warning / Avertissement: Risk of scalding! / Risque de brûlures !**

If you operate your device without a safety hood and at high heating bath temperatures there is a risk of scalding./

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Please wear safety clothes (lab coat, gloves, safety glasses) every time you are in near proximity to the device during operation./

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The heating starts quickly and if empty, the heating bath can be stained./

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Switch on heating only when heating bath is filled sufficiently with fluid./

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## Operation

### Adjust heating bath temperature: Hei-VAP Value, Hei-VAP Digital

- ✓ Device is switched on.
- ✓ Heating bath is filled with water or oil.
- ✓ Evaporating flask is immersed into heating bath.
- Set temperature with temperature knob.
- Press heating key  to start heating.
- ✓ The heating key is illuminated as long as the heating function is on.
- ✓ Control light „heating“ is lit as long as device heats up and set temperature has not been reached yet.
- Press heating key  to stop heating.

### Set heating bath fluid (Hei-VAP Advantage only)

- ✓ Device is switched on.
- ✓ The frame blinks over previously selected medium “oil” or “H<sub>2</sub>O”.
- Select your heating bath fluid within 10 seconds.
- ✓ The frame jumps to the field “temperature”.

After 10 seconds of inactivity, the frame over heating bath fluid disappears.

If you still need to select the heating bath fluid, switch device off and on again.

### Adjust heating bath temperature: Hei-VAP Advantage

- ✓ „bath temp“ is selected.
- Enter desired temperature and confirm.
- Press heating key  to start heating.
- ✓ Heating key is illuminated as long as heating function is on.
- ✓ Actual temperature is shown under „bath temp“.
- Adjust temperature via „bath temp“ field.
- Press heating key  to stop heating.

With heating bath fluid water temperatures above 100°C are not permitted due to safety reasons. Please use oil as heating bath fluid if you adjust temperatures above 100°C.

## ➤ Determine vapor temperature

### (Hei-VAP Advantage only)

- ✓ A vapor temperature sensor is inserted to the condenser and connected with the temperature sensor socket (see chapter "Assembly, Mount device", section "Vapor temperature sensor").
- Switch on device.
- ✓ Actual vapor temperature is shown under „vapor temp“ in the display.

## ➤ Setting duration for process

### (Hei-VAP Advantage only)

Duration of process can be set with the function „timer“.

Adjustable process time: 0 min. – 999 min.

- ✓ Device is switched on.
- ✓ Speed and temperature are set.
- Enter process time in "timer" field and confirm.
- Start heating and rotation.
- ✓ As soon as rotation is started, the timer clock starts to run in reverse.
- ✓ A down arrow is shown under „min“.
- ✓ After end of process time:
  - heating bath turns off
  - frame around time blinks „000“

Arrow pointing down: timer active

### **Variant hand lift**

- rotation is reduced to 20 rpm

### **Variant motor lift**

- rotation switches off
- lift drives up

## ➤ Evaporating and receiving flask

**Warning:**  
**Risk of injury!**

If device is switched on during assembly of evaporating flask, rotation can be started accidentally and might lead to injuries.

Mount and remove evaporating flask only while device is switched off.

**Warning:**  
**Risk of glass breakage!**

If a flask is damaged, glass can break causing serious injuries.

Work with undamaged glass only and check all glassware regularly for any damage.

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
### Mount and remove evaporating flask

#### Mount evaporating flask

- Insert evaporating flask with Easy-Clip open over the vapor tube.
- Tilt Easy-Clip over flange of evaporating flask.
- ✓ Easy-Clip is clicked into place.
- ✓ Evaporating flask is safely assembled.

Easy-Clip open

Easy-Clip closed

#### Remove evaporating flask

Removing the evaporating flask:

- Heating must be switched off.
- Rotation must be at „0“.
- Vacuum pump must be switched off.
- Lift with the evaporating flask must be in top position (see chapter "Start-Up, Heating bath and lift").
- Device must be switched off.

#### **Danger:** **Risk of scalding!**

If the evaporating flask has not cooled down it may cause severe scalding.

- Ensure that the evaporating flask has cooled down prior to exchanging it.
- The evaporating flask should not show a temperature above 50 °C or 122 °F.
- Wear protective clothing.
- Wear appropriate gloves to ensure a safe grip.

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- With 1-2 turns counter clockwise unscrew coupling of the flask clamp Easy-Clip.
- Hold evaporating flask with one hand and open the Easy-Clip.
- Remove evaporating flask carefully.

### Mount and remove receiving flask

#### Mount receiving flask

- Hold receiving flask with open end onto the flange of the condenser.
- Slide clamp over flange of condenser and receiving flask.
- Tighten knurled screw of clamp carefully.
- The receiving flask is mounted.

#### Remove receiving flask

- Hold receiving flask with one hand.
- Loosen knurled screw of clamp.
- Remove clamp from flange.

Knurled screw

Clamp

### ▶ Dispensing sample

Regulate refill of sample during a process through the inlet tube.

- ✓ Glassware is under vacuum.
- Turn inlet valve counter clockwise.
- ✓ Inlet valve is open and sample is drawn into evaporating flask.
- ✓ The further you open the inlet valve, the faster the sample will be drawn in.
- Turn inlet valve clockwise.

Refill of sample will be reduced respectively stopped.

Inlet valve

Ventilation cap

### ▶ Ventilation

The system can be ventilated with the ventilation cap manually.

- To ventilate turn ventilation cap slowly counter clockwise.
- ✓ System is ventilated.
- Turn ventilation cap slowly clockwise.
- ✓ Air inlet is closed gradually.

### > Cleaning

The housing parts and all surfaces of the device allow for cleaning with a damp cloth and if required a mild soap lotion.

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If the device is still connected to power during cleaning, you might get in contact with live parts due to accidentally entering liquids.

Prior to cleaning:

- ensure that the device is switched OFF
- ensure the plug is pulled from the electric socket.

**Caution:**  
**Damage to the device**

Accidentally entering liquids will cause the device to fail .  
Only wipe with a damp cloth. Never use a soaked cloth.

**Caution / Attention:**  
**Damage to the surfaces / Dommages sur les surfaces**

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### Cleaning instructions for heating bath

Over a shorter or longer period of time rust and deposition of solids may be observed.

Rust can be removed easily with stainless steel polish.

Lime deposits can be removed with standard cleaning agents based on acetic or citric acid.

Please follow the manufacturer's instructions in both cases.

### Cleaning instructions for coated glassware sets

High temperatures may damage the coating of your glassware sets.

Do not exceed a maximum temperature of 50° C when washing coated glassware sets.

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Never continue working with a damaged or a worn-off PTFE seal. Otherwise the unit can be damaged accordingly./

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Seals must be checked on a regular basis and must eventually be replaced./

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The device is maintenance-free. Any necessary repair must be performed by an authorized Heidolph distributor.

Please contact Heidolph Instruments or your local Heidolph Instruments distributor.

## Troubleshooting

### ➤ Malfunctions and troubleshooting

General Malfunctions	Possible reason	Troubleshooting
<ul style="list-style-type: none"> <li>Device cannot be switched on.</li> <li>Display doesn't switch on.</li> </ul>	No power supply	<ul style="list-style-type: none"> <li>Check size and shape of plug and compatibility to your electrical socket</li> <li>Check circuit breakers</li> </ul>
	Fuse defective	Replace device fuse (see chapter "Assembly, Electrical Connections")
	ON/OFF-switch off	Switch on ON/OFF-switch
	ON/OFF-switch defective	Please contact your local Heidolph Instruments distributor
No heating function	Connecting cable heating bath not connected	Connect cable
	Heater defective	Please contact your local Heidolph Instruments distributor
	Overheat protection on	<ul style="list-style-type: none"> <li>Heating bath is filled: Cool down heating bath, reset overheat protection (see next section)</li> <li>Heating bath is empty: Contact your local Heidolph Instruments distributor</li> </ul>
No rotation function	Speed controller turned left as far as it will go (Hei-VAP Value, Hei-VAP Value Digital)	Adjust speed
	Drive defective	Please contact your local Heidolph Instruments distributor
Lift not working	Lift is at end of stop	Move lift in opposite direction
	Lift is at height stop	Change height stop
	Mechanic / Motor defective	Please contact your local Heidolph Instruments distributor
No vacuum	Vacuum pump switched off	Switch vacuum pump on
	Vacuum valve defective	Replace vacuum valve
Insufficient vacuum	Leakage in system	<ul style="list-style-type: none"> <li>Check seals and connections</li> <li>Check joints and grease if necessary</li> </ul>
	Vacuum pump defective	Refer to manufacturers instructions
No display of vapor temperature	Vapor temperature sensor not connected	Connect vapor temperature sensor
Device switches off (Hei-VAP Advantage)	Process time activated	Check and turn off timer if necessary

## Troubleshooting

Hei-VAP Value Digital Error message (blinking code)	Possible reason	Troubleshooting
**	Cable connection between controller and base unit disconnected	Connect cable
	Failure in database	Please contact your local Heidolph Instruments distributor
***	<ul style="list-style-type: none"> <li>▪ Heating defective</li> <li>▪ Overheat protection on</li> </ul>	<ul style="list-style-type: none"> <li>▪ Heating bath is filled: Cool down heating bath, reset overheat protection (see next section)</li> <li>▪ Heating bath is empty: Contact your local Heidolph Instruments distributor</li> </ul>
	Safety temperature exceeded by +5°C	Please contact your local Heidolph Instruments distributor
	Heating bath fluid vaporized	<ul style="list-style-type: none"> <li>▪ Switch off device</li> <li>▪ No heating; reset overheat protection (see next section)</li> <li>▪ Refill heating bath fluid</li> </ul>
	Heating bath sensor difference > 10°C	Please contact your local Heidolph Instruments distributor
	<ul style="list-style-type: none"> <li>▪ Temperature sensor heating bath defective</li> <li>▪ Temperature of heating bath fluid &gt; 217 °C</li> <li>▪ Temperature of heating bath fluid &lt; 0 °C</li> </ul>	Please contact your local Heidolph Instruments distributor
****	<ul style="list-style-type: none"> <li>▪ Potentiometer temperature defective</li> <li>▪ Potentiometer speed defective</li> </ul>	Please contact your local Heidolph Instruments distributor
*****	Motor overloaded	<ul style="list-style-type: none"> <li>▪ Switch off device</li> <li>▪ Contact your local Heidolph Instruments distributor</li> </ul>
Display shows --- instead of vapor temperature	Vapor temperature sensor defective	Replace vapor temperature sensor

## Troubleshooting

Hei-VAP Advantage: Error message	Possible reason	Troubleshooting
1	<p><b>Motor lift only:</b> Transportation lock activated</p> <p>Cable connection between controller and base unit disconnected</p>	<ul style="list-style-type: none"> <li>Remove transportation lock</li> <li>Move lift upwards</li> <li>Switch device off and on</li> </ul> <hr/> <ul style="list-style-type: none"> <li>Connect cable</li> <li>Switch device off and on</li> </ul>
3	<ul style="list-style-type: none"> <li>Heating defective</li> <li>Overheat protection on</li> </ul> <hr/> <p>Safety temperature exceeded by +5°C</p> <hr/> <p>Heating bath fluid vaporized</p> <hr/> <p>Heating bath sensor difference &gt; 10°C</p> <hr/> <ul style="list-style-type: none"> <li>Temperature sensor heating bath defective</li> <li>Temperature of heating bath fluid &gt; 217 °C</li> <li>Temperature of heating bath fluid &lt; 0 °C</li> </ul>	<ul style="list-style-type: none"> <li>Heating bath is filled: Cool down heating bath, reset overheat protection (see next section)</li> <li>Heating bath is empty: Contact your local Heidolph Instruments distributor</li> </ul> <hr/> <p>Please contact your local Heidolph Instruments distributor</p> <hr/> <ul style="list-style-type: none"> <li>Switch off device</li> <li>No heating: reset overheat temperature (see next section)</li> <li>Refill heating bath fluid</li> </ul> <hr/> <p>Please contact your local Heidolph Instruments distributor</p> <hr/> <p>Please contact your local Heidolph Instruments distributor</p>
5	Motor overloaded	<ul style="list-style-type: none"> <li>Switch off device</li> <li>Contact your local Heidolph Instruments distributor</li> </ul>
Display shows --- instead of vapor temperature	Vapor temperature sensor defective	Replace vapor temperature sensor

## Overheat protection: reset

- Switch device off. Remove cooled down evaporating flask.
- Disconnect heating bath plug and remove heating bath.
- If necessary empty heating bath and turn it upside down.
- Push the red knob on the bottom of the heating bath with a pointed object.
- ✓ Thermal cut-out is reset.

Reset overheat protection

If you experience a malfunction which can not be resolved, please contact your authorized Heidolph distributor immediately.

## Electrical connections

- Fuses must **only** be replaced by a **professional electrician**.
- **Repairs** of any kind are allowed by **qualified professional electricians only**. Any improper repair can result in a dangerous situation. Contact your local Heidolph Instruments distributor for any repair you may have.

### Warning: Electric shock!

This device must only be connected to a grounded electrical socket.



The unit and voltage must match. The rating plate on the back side of the unit provides all voltage specifics (see picture below).

When shipped, the unit is grounded. In case the original plug is changed, the new plug must have a protective conductor!

If the device is still connected to power while changing the fuse, you might get in contact with live parts.

Prior to replacing the fuse ensure that the unit is switched off and pull the plug from the electric socket. Only use original fuses from Heidolph.

After changing a fuse check the device for safe operating condition according to IEC 61010-1.

### Color code for electrical connections:

Europe	
<b>GREEN/YEL-LOW</b>	PE: Protective conductor ( <b>Earth</b> )
<b>BLUE</b>	N: Neutral conductor
<b>BROWN</b>	P: Phase

North America	
<b>GREEN</b>	PE: Protective conductor ( <b>Earth</b> )
<b>WHITE</b>	N: Neutral conductor
<b>BLACK</b>	P: Phase

The unit is connected with the power cord supplied to the electrical socket. It plugs into the housing on the rear side of the unit.

For countries utilizing a plug other than the standard one supplied:

- The plug supplied may only be changed by a professional electrician.
- If utilizing an adaptor ensure that it is grounded and approved by the local regulations.
- Wiring, adaptor and plug must at least meet the electrical power stated on the rating plate.

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


## Transportation lock

(Motor lift only)

Transportation lock must be removed before you start. Keep the transportation lock for possible return shipments.

The transportation lock consists of three M5x8 screws and a connection plate with drilled holes to fix the position.

- ✓ Base unit is placed on a stable surface.
- Press and hold lift contact on the lift.

Lift contact

Transportation lock

- Loosen the three screws with allen key supplied and remove plate.
- Release lift contact.
- ✓ Lift drives up automatically.
- Keep all individual parts for subsequent use.


Transportation lock

## Heating bath

### Attach heating bath

- ✓ A movable base plate with an enhanced guide rail is located on the right side of the base unit.
- Place heating bath with gap on the guide rail.
- Move heating bath carefully to the left and right to lock it.

### Move heating bath

It may be necessary to increase the distance between heating bath and drive when bigger evaporating flasks or spacers between evaporating flask and vapor tube are used.

- ✓ Base unit is placed on a stable surface.
- Allow sufficient space to move heating bath.
- ✓ On the right side the heating bath slides out approximately 20 cm to allow for adjustment.
- Hold heating bath on both handles and slide it to the right until you have reached the desired position.

## Vapor tube



### Assembly parts

- Vapor tube
- Flask clamp Easy-Clip
- Coupling vapor tube
- Clamping sleeve

The flask clamp Easy-Clip is pre-assembled in NS 29. For smaller flasks the Easy-Clip NS 24 is enclosed.

The device is delivered with pre-assembled coupling for the vapor tube and the condenser.


- Hold the locking knob beneath "lock" for a quarter turn.
- Unscrew coupling of the vapor tube along with the vapor tube on the right side of the drive.
- Unscrew coupling of the condenser on the left side of the drive.

## Assembly

- Remove tension spring and PTFE seal on the condenser side.

Tension spring

PTFE-seal

**Caution:**  
**Possible contamination, leakage, damage to your unit**

Due to a damaged PTFE seal your SET vacuum may not be reached, your sample may be contaminated, your unit may be damaged.

Especially the sealing lip of the vacuum seal is susceptible to damage in handling. Never operate your system with a damaged or worn-off PTFE seal! Always check the condition of your PTFE seal regularly and if required replace it.

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- ✓ Coupling and Easy-Clip are screwed together.
- Push coupling of the vapor tube along with Easy-Clip on the vapor tube (direction and order as shown under "Assembly parts").
- Push clamping sleeve with the bead to the glass on vapor tube until it audibly locks.
- Insert vapor tube with clamping sleeve first from below to the drive and screw it tight.
- At the very last quarter turn, press locking knob.
- ✓ Coupling locks in and vapor tube is securely mounted.

Vapor tube, clamping sleeve and coupling

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Work with undamaged glass only and check all glassware regularly for any damage./

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Heidolphs coated glassware sets aim to minimize the risk of injury due to glass breakage. However they are only suitable for temperatures up to 80° C.

### **PTFE seal**

- Slide PTFE seal with the labeling "Motor Side" facing drive onto the vapor tube.
- ✓ The PTFE seal is tightened as far as it will go.

PTFE seal

### **Condenser: Glassware set G1 and G3**

- Push coupling of glassware over condenser flange.
- Push tension ring over condenser flange.
- Screw condenser to the drive hand-tight.

Tension ring      Coupling

- ✓ The condenser is mounted to the drive.

### Condenser:

#### Glassware set G5 dry ice

- ✓ The condenser is mounted to the drive.
- Screw inlet tube with inlet valve onto opening opposite of vapor tube (see chapter "Assembly, Inlet tubes").
- If using inert gas screw inlet tube with inert gas valve onto opening above inlet valve.

Inert gas valve

Inlet valve

- Place red seal ring into notch at top of condenser flange.
- Place white PTFE ring with collar up on to condenser flange.
- Insert the jacket trap into condenser.
- ✓ Due to vacuum inside condenser, the flange is drawn in and the system is tight.
- To protect the dry ice, place black cap on top.

## Assembly

### Condenser: Glassware set G6 (Condenser for reflux distillation)

Additional assembly parts on the vapor tube:

- PTFE seal ring
- Spacer sleeve
- Additional coupling

✓ Vapor tube together with PTFE seal is mounted.

→ Insert PTFE seal ring into spacer sleeve.


→ Place spacer sleeve including PTFE seal ring with the narrow side facing to vapor tube on condenser side.

→ Screw on additional coupling with external thread pointing left.



→ Tighten distributor with coupling including tension spring on external thread.

→ Insert condenser to the upper joint of distributor.

→ Screw inlet tube with inlet valve onto opening opposite vapor tube (see chapter "Assembly, Inlet tubes").

Distributor

Distributor

Inlet tube


## Assembly

- ✓ Distributor is mounted.
- Turn reflux valve clockwise until O-ring is clearly visible and white plug blocks output to receiving flask.
- ✓ Valve is closed and sample flows back to evaporating flask.
- Turn reflux valve counter clockwise until white plug unblocks the output of receiving flask.
- ✓ Valve is open and sample flows into receiving flask.

(For inlet valve and inlet tube see chapter "Assembly, Inlet tubes".)

The system is not vacuum-tight if valve is opened too far. Open valve only in a way that the white seal ring beneath knurled screw remains visible 0.5 cm wide.


## Assembly

### Condenser support for vertical condenser: Glassware set G3 and G6

#### Assembly parts:

- Support rod
- Boss head with condenser clamp
- Allen key

Knurled screw condenser clamp    Boss head

Support rod    Condenser clamp    Allen key

On the left side of the unit a drilling hole is placed to fix the support rod.

- Insert support rod with the flattened side downwards into hole.
- Secure support rod with allen screw using supplied allen key.
- Slide boss head of cross pole over support rod.
- Fix cross pole with knurled screw at support rod.
- Fix the condenser clamp with knurled screw around the condenser.

Hole

Allen screw

Knurled screw boss head

### Condenser support for vertical condenser: Glassware set G5

#### Assembly parts:

- Support rod
- Boss head with tension band
- Allen key
- Mount condenser support analog to condenser support for G3 and G6.
- Fix condenser with the tension band.

Tension band    Allen key

Knurled screw    Boss head

## Inlet tubes


### Assembly parts:

- 2x Inlet tube
- 1x Inlet valve

On the bottom left side of the condenser a ground opening is located to include stopcock with inlet tube.


### Inlet tube 1

- Shorten supplied PTFE tube so that the end of the tube leads to the evaporating flask and attach it to the stopcock.
- Slide drain disk onto PTFE tube so that it is placed below the vapor opening of condenser.
- ✓ Condensed fluid cannot flow back into evaporating flask through inlet tube.
- Unscrew screw cap from condenser.
- Slide seal ring and screw cap approx. 3 cm onto PTFE tube.
- Insert PTFE tube with seal ahead into inlet opening and tighten screw cap.


### Inlet tube 2

- Unscrew the screw cap including the seal ring of inlet opening.
- Slide seal ring and screw cap approx. 3 cm onto PTFE tube.
- Insert PTFE tube with seal ahead into inlet opening and tighten screw cap.
- ✓ Inlet tubes are mounted.

### Inlet tubes glassware set G1

Openings for inlet tubes are located on the top of the condenser.

## Vapor temperature sensor, AUTOaccurate sensor

### Vapor temperature sensor

(Optional accessory, Hei-VAP Advantage only)

### AUTOaccurate sensor

(Glassware set G3 and G6 and upgrade to Hei-VAP Precision only)

Opening for vapor temperature sensor and AUTOaccurate sensor is located at the bottom of the condenser.

Opening vapor temperature sensor

- Remove sealing cap of condenser.
- Insert vapor temperature sensor into condenser.
- Tighten vapor temperature sensor with screw cap.
- Plug in connector plug with white mark upwards into receptacle on the left side of unit.

Connector plug vapor temperature sensor


Screw cap vapor temperature sensor

## Evaporating flask angle and depth

### Evaporating flask angle

Evaporating flask angle can be adjusted by turning the drive.

- ✓ Rotation is set to 0 rpm.
  - ✓ Evaporating flask is mounted (see chapter "Operation", section "Mount evaporating flask").
- Hold condenser with your left hand.
  - Slide locking mechanism on lift pillar to the right and hold it against the spring force.
  - Move condenser carefully until connected evaporating flask has the correct angle.
  - Release locking mechanism.
  - If necessary move condenser slightly until it locks in.

Locking

Lift pillar

### Flask depth

A stop bolt is located on the right side of the lift pillar with which you can adjust a constant depth for the evaporating flask.

- ✓ Rotation is set to 0 rpm.
- ✓ Evaporating flask is mounted (see chapter "Operation", section "Mount evaporating flask").
- ✓ Height stop is at top position.

#### **Adjust height without guard hood:**

- Move down lift until it has reached the desired depth (see chapter "Start-Up, Heating bath and lift").
- Turn fixing screw on stop bolt of lift pillar approx. 180° counter clockwise.
- Slide stop bolt via fixing screw downwards as far as it will go.
- Fasten fixing screw clockwise.
- Move the lift upwards.
- ✓ The fixing screw allows to individually adjust the lowest flask position.
- ✓ The lift cannot be lowered beyond this point.

#### **Adjust height with guard hood:**

- Hold guard hood with one hand.
- Move guard hood manually when adjusting the lift.

Fixing screw

Lift contact

Stop bolt

## Control panel

Control panel can be removed from console to be fixed to a wall or fume hood.  
Connection cable is 1 meter long.

### **Warning: Risk of poisoning!**

If operating with control panel inside of fume hood, there is a risk of harm if fume hood is opened during process.

Remove control panel and place it outside the fume hood.

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### **Open cable carrier**

- Grab control panel on the top left side and tip it to the front.
- Avoid pressure to the LCD display.
- ✓ The cable winding is surrounded by a shell.
- Pull shell of cable winding open.
- Wrap cable to the desired length on or off.

### **Position control panel**

- If necessary turn shell 180°.
- Press shell and control panel together until they lock audibly.
- Place control panel on a clean and stable surface or attach it to a wall.

### **Insert control panel into console**

- Turn shell into initial position with the wide side downwards.
- Wrap cable up to the yellow colored cable mark.
- Press shell and control panel together.
- If necessary plug in connecting plug with white colored tag upwards into console.

Yellow cable mark

- Place control panel with an angle of approx. 30° forward into the console.
- Press control panel into console until it locks audibly.

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Flasks which show even the slightest signs of damage may burst and depending on the solvent poison the environment./

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Increase vacuum pressure and/or leakage pressure only if necessary. Do not increase pressure/under-pressure beyond the limit required by your solvents./

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To connect the vacuum you need:

- a vacuum source
- a vacuum tube (Ø 7-8 mm inside)

#### Glassware set G1, G3 and G6

The opening of the vacuum connection is located on the top of condenser

- Unscrew either one of the screw caps with nipple.
- Push vacuum tubing onto nipple.
- Tighten screw cap with nipple and tubing.
- Connect other end of vacuum tube with vacuum source.

#### Glassware set G5

The opening for vacuum connection is located on the top of condenser.

Vacuum connection G1, G3 and G6

Screw caps vacuum

Vacuum connection G5


Opening vacuum connection

## Assembly

### Vacuum connecting scheme (demonstration)

## Connect cooling water

**Warning:**  
**Risk of explosion**

A closed valve in the recirculation system will build up high pressure in the glassware and may result in glass bursting.


The pressure in the condenser shall not exceed 2 bar to ensure that a higher pressure never occurs.

No valve or restriction may be built into the recirculation system.

**Assembly parts:**

- 2x tube nipples  
(included in scope of delivery)
- 1x tube, inner Ø 7-8mm  
(not included in scope of delivery)


Tube nipples

Input and output for cooling water is located on the back of condenser.

- Remove screw cap from input and output of the cooling water.
- Push screw cap over tube nipples so that the black seal is placed inside the cap.
- Push intake and outflow hoses for cooling water through screw cap over one of each tube nipple.
- Tighten screw cap with seal side up to the openings of condenser.
- ✓ Intake and outflow of cooling water are connected.

Input and output cooling water

## Mount tubes

To mount tubes an optional tubing clamp is available.


- Install tubing clamp on the back of device.

## ➤ Additional accessories


### Guard hood

- Fix carrier guard hood at height stop with two allen screws supplied.
- Release clamping screw and align guard hood while moving it horizontally.
- Tighten clamping screw
- ✓ Guard hood is mounted.
- Open guard hood with handle.

### Safety shield

- Move lift upwards.
- Pull up safety shield to heating bath and place it so that the guard hood can be opened easily, where applicable.
- Turn eccentric screw 180° and fix safety shield to heating bath.

## ➤ Dismantling, transportation and storage

### Dismantling

**Warning:**  
**Risk of scalding!**

Hot glass surfaces and heating bath fluids may result in scalding prior to dismantling.

Ensure that all parts have cooled down to room temperature.

**Risk of poisoning!**

Minor leakages may result in solvents and toxic substances being released into the environment.

Ensure that fluids and solvents are collected in appropriate vessels.

Ensure that the surface of the unit is not contaminated with solvents and hazardous materials.

To clean the unit you may run a distillation with harmless substance.

**Risk of electric shock!**

If the device is still connected to power during dismantling, you might get in contact with live parts.

Prior to dismantling unit ensure that the unit is switched off and the plug is pulled from the electric socket.

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- ✓ The unit is switched off.
- Remove all connections to other devices.
- Remove all glass ware and heating bath.
  - Motor lift**
  - Switch device ON.
  - Move lift downwards as far as it will go.
  - Screw transportation lock while holding down lift contact.
  - Switch the device OFF.
- Disconnect plug from power socket and device socket.

## Transportation and storage

- Store and transport the unit and its components only if they were emptied, cleaned and dismantled.
- To carry the device wear safety gloves and lift the unit from below.
- Store and transport the device and parts of device in the original packing material or alternatively in an adequate container to prevent damage.
- Seal the package carefully against unauthorized or accidental opening.
- Store the device in a dry and frost-free place.

Improper storage and transportation may result in damages to the system and the mechanical components.

Avoid any kind of shocks during transportation.

## ➤ Scope of delivery

Component	Variant	Quantity	Product number
Model with hand lift and glassware set G3: Hei-VAP	Value	1	560-01300-00
	Value Digital	1	560-01302-00
	Advantage	1	561-01300-00
Model with motor lift and glassware set G3: Hei-VAP	Advantage	1	562-01300-00
PTFE/FKM seal		1	23-30-01-01-30
Vapor tube		1	514-00000-01
Clamping sleeve		1	23-30-01-05-31
Transportation lock (motor lift)		1	11-300-006-28
Coupling		1	23-09-03-01-03
Tension ring		1	22-03-02-01-05
Flask clamp Easy-Clip NS 29		1	23-30-01-05-29
Flask clamp Easy-Clip NS 24		1	23-30-01-05-57
Joint clamp		1	515-42000-00
Operating Manual English / German		1	01-005-004-79
Warranty registration / Confirmation of condition		1	01-006-002-78
EU Declaration of Conformity		1	01-001-025-01
Power cord		1	depending on country

## ➤ Accessories

Component	Quantity	Product number
Device fuses (115 V and 230 V)	2	11-300-009-40
Vapor temperature sensor	1	569-00030-00
Guard hood	1	569-00010-00
Tube (suitable for vacuum and water)	1	591-35000-00
PTFE/FFKM seal	1	23-30-01-06-70
Glassware set		see general catalog
Vacuum pumps and chillers type and size vary		see general catalog

Please observe the following **permissible temperature range** when using **coated glassware sets**:

- G5 (dry Ice glassware set): -76° C bis 150° C
- any other glassware sets: -40° C bis 80° C

Additional accessories are shown in the general catalog or on our website at [www.heidolph.com](http://www.heidolph.com)

## ➤ Technical data

### Hei-VAP Value, Hei-VAP Value Digital, Hei-VAP Advantage

Standard supply voltage	230 V/50 Hz or 115 V/60 Hz	
Supply power (W)	1400	
Protection class (DIN EN 60529)	<ul style="list-style-type: none"> <li>▪ IP 20 device</li> <li>▪ IP 67 cable heating bath</li> </ul>	
Lift distance (mm)	155	
Drive	EC motor	
Rotation speed (rpm)	10 - 280	
Speed control	electronic	
Heating bath	<ul style="list-style-type: none"> <li>▪ inner Ø: 253 mm, outer Ø: 291 mm</li> <li>▪ Material: V4A-steel (1.4404) stainless steel X2CrNiMo 17-12-2</li> <li>▪ Volume 4,5 L</li> </ul>	
Heating capacity (W)	1,300	
Temperature range (°C)	20 - 100 H <sub>2</sub> O / 20 - 210 oil	
Temperature control	electronic	
Temperature accuracy Heating bath (°C)	±1	
Overtemperature protection heating bath	cut off at 5 °C difference to set temperature	
Overheat protection (°C)	250	
Condensing surface (cm <sup>2</sup> )	1,400	
Airborne noise level	< 85 db(A)	
Evaporation rates (L/h) ΔT* = 40°C (continuous operation)	<ul style="list-style-type: none"> <li>▪ Toluene 8.5</li> <li>▪ Acetone 5.8</li> <li>▪ Ethanol 3.5</li> <li>▪ Water 1.2</li> </ul>	
<b>Lift</b>	<b>Manual</b>	<b>Motor</b>
Weight (without glassware) (kg)	14.7	15
Lifting speed (mm/s)	-	25
Dimensions with glassware set G3 (w x d x h) (mm)	739 x 490 x 887	739 x H 420 x 887
<b>Device</b>	<b>Value Digital</b>	<b>Advantage</b>
Display bath temperature (°C)	digital	digital
Display vapor temperature (°C)	-	digital (required: vapor temperature sensor)
Display rotation speed	-	digital
Timer	-	yes

\* ΔT = Difference between heating bath temperature and boiling temperature

## ➤ Contact/Technical Support

### Questions / Repair work

If any aspect of installation, operation or maintenance remains unanswered in the present manual, please contact us.

For repairs please call Heidolph Instruments or your local authorized Heidolph distributor.

---

#### **Warning: Danger of poisoning!**

Contaminated devices can lead to severe injury or death of our employees!

When shipping items for repair that may have been contaminated by hazardous substances, please

- advise exact substance
  - take adequate protective action towards our parts receiving and service personnel
  - mark the pack in accordance with Ordinance on Hazardous Substances
- 

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Included at the end of this manual you will find a „Confirmation of condition“.

→ Kindly copy and fill in this form and submit it prior to shipping the device for repair.

## Service

### Contact details:

---

#### **Heidolph Instruments Germany**

Heidolph Instruments GmbH & Co. KG  
Technical Service  
Walpersdorfer Str. 12  
D-91126 Schwabach / Germany

Tel.: +49 - 9122 - 9920-74  
Fax: +49 - 9122 - 9920-84

E-Mail: [service@heidolph.de](mailto:service@heidolph.de)  
[www.heidolph.de](http://www.heidolph.de)

#### **Heidolph Instruments North America**

Phone: 1-866-650-9604  
E-mail: [service@heidolph.com](mailto:service@heidolph.com)  
[www.heidolphNA.com](http://www.heidolphNA.com)

#### **Heidolph Instruments United Kingdom**

Phone: 01799 - 5133-20  
E-mail: [service@radleys.co.uk](mailto:service@radleys.co.uk)  
[www.heidolph-instruments.co.uk](http://www.heidolph-instruments.co.uk)

#### **All other countries**

You will find contact details of your local Heidolph distributor at  
[www.heidolph.com](http://www.heidolph.com)

---

## ➤ Warranty

---

Heidolph Instruments provides a 3-year warranty for the products described here (excluding glass and wearing parts) if you register using the warranty card enclosed or online ([www.heidolph.com](http://www.heidolph.com)). The warranty is valid from the point of registration. The serial number is also valid without registering.

The warranty covers part and manufacturing defects.

---

In the case of a part or manufacturing defect, the device shall be repaired or replaced free of charge under the terms of the warranty.

Heidolph Instruments shall not assume liability for any damages incurred as a result of improper handling or transport.

#### Warranty claim?

- ➔ Please inform Heidolph Instruments or your local Heidolph distributor should you wish to make a warranty claim.
-

## ➤ Confirmation of condition

→ In the case of repair, copy and complete the Confirmation of condition and send it to Heidolph Instruments.

### 1. Details about the device

Product number \_\_\_\_\_

Serial number \_\_\_\_\_

Reason for repair \_\_\_\_\_

\_\_\_\_\_

### 2. Has the device been cleaned, decontaminated/sterilized?

Yes \_\_\_\_\_ No \_\_\_\_\_

### 3. Is the device in a condition which does not represent any health threats for the staff of our service department?

Yes \_\_\_\_\_ No \_\_\_\_\_

If not, which substances has the device come into contact with?

\_\_\_\_\_

### 4. Legally binding declaration

The customer is aware of being legally liable to Heidolph Instruments for any damages arising from incomplete and incorrect information.

Date \_\_\_\_\_

Signature \_\_\_\_\_

### Company stamp

#### Please note

\_\_\_\_\_

The shipper is responsible for the return of the goods in well packed condition, suitable for the mode of transport.

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This is a translation of the original operating manual. We reserve the right to make any technical modifications without prior notification. This documentation is not subject to revision; the latest version of this documentation can be found on our homepage. The copyright for all texts and images is held by Heidolph Instruments GmbH & Co. KG.
© Heidolph Instruments GmbH & Co. KG