Automatic module

Hei-VOLUME Distimatic Pro





- ✓ Saves time
- ✓ Increases throughput
 - Reduces costs



Saves time. Increases throughput. Reduces costs.

Dealing with large volumes does not have to be time-consuming: The Hei-VOLUME Distimatic module allows you to fully automate the filling and emptying of your Hei-VAP Industrial large-scale rotary evaporator.

Reduced work load through automation

- Automatic filling and refilling of the rotary flask
- Automatic distillate drainage through overpressure
- Automatic overpressure drainage of the residue (the 24/7 models)
- Automatic shutdown of all system components at the end of the process and in the event of a disturbance
- Automatic cleaning program can be selected











Very high usability

- Intuitive operation with push rotary knobs and 7" touchscreen
- Guided menu navigation for the smooth setting and storage of process parameters
- Process flow display with pictograms and colors for clear identification of process progress
- Removable control panel for process adjustment even from remote locations
- Multilingual user interface

Complete efficiency and safety through smart sensor technology

- Increased distillation rates compared to manual operation by maintaining constant conditions
- Calibration-free sensors reliably detect any medium
- Individual process setting through selectable operating modes (sensor, time, manual)
- Safe operation day/night even unattended
- Flexible configuration with accessories for even more safety
- Space-saving solution for high-volume processing
- Cost-effective alternative to purchasing larger evaporator systems







Can be combined with all Heidolph large-scale rotary evaporators



Available in Hei-VAP³ packages – by a factor 3 more efficient evaporator system compared to manual operation

Hei-VOLUME Distimatic Pro Industrial for glassware A 591-61210-00

Hei-VOLUME Distimatic Pro 24/7 Industrial for glassware A 591-61110-00

Hei-VOLUME Distimatic Pro Industrial for glassware R 591-61220-00

Hei-VOLUME Distimatic Pro 24/7 Industrial for glassware R

Hei-VAP³ Pro Industrial High Boiling Package

Hei-VAP³ Pro Industrial Solvent Recovery Package 518-14162-00

Hei-VAP³ Pro Industrial 24/7 High Boiling Package 518-11165-00

Hei-VAP³ Pro Industrial 24/7 Solvent Recovery Package 518-14164-00



