

# Precise

## Dosing and Dispensing

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.

Safety for continuous operation: the motor will be switched off if a high thermal load situation occurs to **prevent accidents**



MADE IN  
**GERMANY**

All models comply with **protection class IP 55** to prevent corrosion and short circuits



**High precision** even for flow rates starting at 0.005 ml/min

Optional foot-pedal remote control via cable to operate pump drive from **outside closed fume hoods** facilitate operation

All models feature a smooth start operation which **prevents spills and splashing media**. The speed will ramp up slowly until your set rpm has been reached

For the use in cell biology just use a pump head which features convex rollers and **does not squeeze** the tubing as conventional rollers