

PRECISE HEATING & STIRRING HEI-PLATE MIX 'N' HEAT MAGNETIC STIRRER

HIGH-RESISTANT HEATING PLATE

REFRACTORY HOUSING

800 W

MADE IN
GERMANY

3-year warranty on
all devices and an op-
erational lifespan of
more than 10 years

Hei-PLATE Mix 'n' Heat Expert

Interface for more precision and comfort

For the Pt1000 or the double sensor Pt1000 to measure the temperature exactly until inside the liquid

USB B interface for a precise and quick reproducibility of your results via our free Hei-Control software

Well-designed features for a carefree working

Different heating modes (FAST, PRECISE, PRECISE+) and safety features handle your samples optimally

Three-stage setting of the stirring speed (start-up behavior). Additionally options for the rotational direction and the stir bar detector

Intuitive handling for a quick start

Standby button in view to switch the device easily on and off

Display per main function (stirring and heating) to have an eye on all values and settings at anytime

The locking function of the push-and-rotary controls protects against accidental changing of the values

LED rings for the easy identification of the active functions and residual heat display

Timer per main function can be set independently
Shows the remaining process duration to the nearest second.

User-friendly menu navigation for the activation/deactivation of well-designed features
Language versions: German, English



ACCESSORIES RECOMMENDATION



Silicone protective cover



Double temperature sensor



Heating and bath attachments



High-performance condenser