

Intelligent Evaporation



MADE IN
GERMANY

The average operational **lifespan of 10 years** is backed by a **3 year warranty** and makes your purchase a worthwhile investment.



Ground-free condenser:
Ventilation cap with
PTFE inlay and threaded
connection **prevents
stuck glass**

Up to 40 % higher evaporation rate with a cooling
surface of 2,200 cm² for all XL glass condensers

Highly resistant and durable PTFE vacuum seals
for many years of continuous operation

The right operating concept
for all requirements:
**with digital display or
innovative touch display.**
Also available as Control
model with Control-Box
for additional control of
vacuum and Hei-CHILL
recirculation chiller

Patented connection system
with Easy-Clip to gently
screw-off flasks and **vapor
tube with clamping sleeve**
prevent stuck glass parts

**The 210 °C safety heating
bath** with ergonomic handles
and pour spout shuts down
when running dry or at
over temperature > 5 °C

For quick access, all models
also have separate knobs
with **LED rings showing
the active functions**

Cable connections with
protection class IP 67
protect against corrosion
and short-circuits

For unlimited evaporation: In conjunction with
the Distimatic automatic module, you can distill
indefinitely and automatically in 24 / 7 operation. The
affordable alternative to large rotary evaporators.

With standby function and
residual heat warning at heating
bath temperatures above 50 °C

Metal connection to the
device console **prevents
tilting of the bath**