

Reax 20 - for 4 bottles (1-16 rpm)

P/N: 541-20004-00



The resilient one with a capacity of up to 4 bottles

The Reax 20 is an ideal shaker for processing and storing samples that tend to sediment. The fields of application range from environmental technology to paint processing and quality assurance in the materials industry.

The maximum permissible weight is 30 kg, so that samples with a high solid content can also be properly mixed. Since this device sets great forces in motion, it is equipped with an emergency stop switch for your safety.

This overhead shaker can be ordered with a capacity of [4](#), [8](#) or [12](#) bottles. Also suitable for mixing cylinders or wide-neck bottles.

Reax 20 - for 4 bottles (1-16 rpm) - Technical Data

| | |
|--------------------------------|--|
| Overheat protection | self-resetting |
| Permissible ambient conditions | 5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity |
| Weight | 23 kg |
| Protection class IEC 60529 | IP 21 |
| Timer | - |
| Rotation speed range | 1 - 16 rpm |
| Power input | 280 W |
| Dimensions (w/d/h) | 490 x 520 x 465 mm |
| Max. Load | 30 kg |
| Operating Mode | - |
| Rotation speed setting | electronic control |
| Orbit / Stroke | - |
| Platform size | - |
| Accessories included | Tension plate for cap dia. up to 77 mm |
| Angle | - |
| Motion | overhead |