

## Silicone two-stop tubing for Cassette small - wt: 0.9 | id: 2.79 | od: 4.59

P/N: 525-30016-00



### Two-stop tubing (40 cm) for Cassette small

- Application: Platinum-coated silicone hose for use in pharmaceuticals and biology
- Features:
  - Extremely smooth interior prevents bacterial growth
  - Biocompatible, minimal adsorption and absorption
  - Best flow properties, high temperature stability
  - Absolutely inert, softener-free
- Material: Polydimethylsiloxane with siliceous earth and silicone additives; translucent white; excellent resistance to initial pressure
- Complies with the standards: USP Class VI meets GLP and NSF
- Temperature range: -80 to +200 °C
- Sterilization: Can be autoclaved for 30 min at 1 bar or sterilized by irradiation
- Restriction: Not suitable for concentrated solvents, oils, acids or dilute caustic soda; relatively high permeability to gas

Used with:

Acids:	conditional
Lyes:	conditional
Solvents:	unsuitable
Pressure:	satisfactory
Vacuum:	good
Viscous media:	satisfactory
Sterile media:	excellent

## Silicone two-stop tubing for Cassette small - wt: 0.9 | id: 2.79 | od: 4.59 - Technical Data

Inner diameter	2.79 mm
Outer diameter	4.59 mm
Wall thickness (wt)	0.9 mm
Max. Pressure (continuous / short time)	0.5 / 1.5 bar
Suction height	7 mH <sub>2</sub> O
SP quick Advantage 06 / Precision 06 (min/max)	not suitable
SP quick Value 06 (min/max)	not suitable
SP quick Advantage 01 / Precision 01 (min/max)	not suitable
SP quick Value 01 (min/max)	not suitable
SP standard/vario Advantage 06 / Precision 06 (min/max)	not suitable
SP standard/vario Value 06 (min/max)	not suitable
SP standard/vario Advantage 01 / Precision 01 (min/max)	not suitable
SP standard/vario Value 01 (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump head C 12 - Cassette small (min/max)	0.29 – 9 ml/min
Hei-FLOW Advantage 01 / Precision 01 - pump head C 4 - Cassette small (min/max)	2 – 36 ml/min
Hei-FLOW Value 01 - pump head C 12 - Cassette small (min/max)	0.69 – 8 ml/min
Hei-FLOW Value 01 - pump head C 4 - Cassette small (min/max)	4 – 37 ml/min
Hei-FLOW Value 01 - pump head C 8 - Cassette medium (min/max)	not suitable
Hei-FLOW Value 01 - pump head C 8 - Cassette large (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable

head C 8 - Cassette medium (min/max)

Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable
head C 8 - Cassette large (min/max)	

Flow rates pertain to	Tygon® (Standard) tubing and water
-----------------------	------------------------------------