

Hei-Vac Valve Industrial

P/N: 591-07210-00



Includes secondary condenser and air intake separation vessel

This vacuum system can be applied in many different fields in order to evacuate, evaporate and pump out gases and vapors.

Thanks to the emission condenser it is possible to achieve a chemical recovery of almost 100%

- Excellent compatibility with chemicals and condensate
- Excellent ultimate vacuum of 2 mbar
- Very quiet and very low vibration
- Very high suction capacity of 3.0 m³/h

Hei-Vac Valve Industrial - Technical Data

Weight	19.9 kg
Diaphragm pump	yes
Suction capacity	3 m³/h
Ultimate vacuum	2 mbar
Combinable with a condenser	yes
Power input	250
Dimensions without condenser (w/d/h)	350 x 275 x 495 mm