

## Viton® two-stop tubing for Cassette small - wt: 0.9 | id: 1.42 | od: 3.22

P/N: 525-50015-00



### Two-stop tubing (40 cm) for Cassette small

- Application: Excellent acid resistance at high temperatures
- Features:
  - Low gas permeability
  - Resistant to solvents and corrosives at high temperatures
- Material: Fluorocarbon rubber, thermoformed Viton® B (67% fluorinated); opaque black
- Complies with the standards: GLP conform
- Temperature range: -30 to +205 °C
- Sterilization: Sterilization is not recommended
- Restriction: Limited performance life

Used with:

Acids:	excellent
Lyes:	excellent

## Premium Laboratory Equipment

Solvents:	varies; tests recommended
Pressure:	good
Vacuum:	good
Viscous media:	good
Sterile media:	satisfactory

## Viton® two-stop tubing for Cassette small - wt: 0.9 | id: 1.42 | od: 3.22 - Technical Data

Inner diameter	1.42 mm
Outer diameter	3.22 mm
Wall thickness (wt)	0.9 mm
Max. Pressure (continuous / short time)	0.5 / 1.5 bar
Suction height	7 mH <sub>2</sub> O
SP quick Advantage 06 / Precision 06 (min/max)	not suitable
SP quick Value 06 (min/max)	not suitable
SP quick Advantage 01 / Precision 01 (min/max)	not suitable
SP quick Value 01 (min/max)	not suitable
SP standard/vario Advantage 06 / Precision 06 (min/max)	not suitable
SP standard/vario Value 06 (min/max)	not suitable
SP standard/vario Advantage 01 / Precision 01 (min/max)	not suitable
SP standard/vario Value 01 (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump head C 12 - Cassette small (min/max)	0.10 – 3 ml/min
Hei-FLOW Advantage 01 / Precision 01 - pump head C 4 - Cassette small (min/max)	0.60 – 14 ml/min
Hei-FLOW Value 01 - pump head C 12 - Cassette small (min/max)	0.23 – 3 ml/min
Hei-FLOW Value 01 - pump head C 4 - Cassette small (min/max)	1 – 15 ml/min
Hei-FLOW Value 01 - pump head C 8 - Cassette medium (min/max)	not suitable
Hei-FLOW Value 01 - pump head C 8 - Cassette large (min/max)	not suitable
Hei-FLOW Advantage 01 / Precision 01 - pump	not suitable

head C 8 - Cassette medium (min/max)

Hei-FLOW Advantage 01 / Precision 01 - pump not suitable  
head C 8 - Cassette large (min/max)

Flow rates pertain to

Tygon® (Standard) tubing and water