

## Hei-VAC Valve Control

P/N: 591-00130-00



- Suction capacity for up to three rotary evaporators at the same time
- Vacuum can be controlled via valve operated or manual vacuum controllers
- Depending on your application you can switch on and off the vacuum pump via switchbox
- A Vacuum valve is necessary, if you want to combine this vacuum pump with the Hei-VAP Expert Control / Hei-VAP Ultimate Control rotary evaporators

Condenser is not in scope of delivery.

# Hei-VAC Valve Control - Technical Data

Weight	12.8 kg
Diaphragm pump	Two-stage
Suction capacity	2 m³/h
Ultimate vacuum	7 mbar
Combinable with a condenser	yes
Power input	180
Dimensions without condenser (w/d/h)	195 x 245 x 310 mm