

## Hei-VAC Vario Station

P/N: 591-00142-00



Fully controllable stand-alone pumping unit including vacuum controller.  
For Hei-VAP Value and Hei-VAP Advantage

- Precise vacuum control avoids bumping and eliminates foaming of your evaporation solution
- Due to automatic vacuum supply to the process parameters you have more time for important laboratory activities

Condenser is not in scope of delivery.

## Hei-VAC Vario Station - Technical Data

Weight	6 kg
Diaphragm pump	Three-stage
Suction capacity	1.7 m <sup>3</sup> /h
Ultimate vacuum	5 mbar
Combinable with a condenser	yes
Power input	160
Dimensions without condenser (w/d/h)	294 x 192 x 299 mm