

Hei-VAC Vario Tec

P/N: 591-00171-00



For Hei-VAP Expert Control / Hei-VAP Ultimate Control

- Recommended for solvents with low or medium boiling points

Condenser is not in scope of delivery.

Hei-VAC Vario Tec - Technical Data

Weight	4.3 kg
Diaphragm pump	Two-stage
Suction capacity	1 m ³ /h
Ultimate vacuum	12 mbar
Combinable with a condenser	yes
Power input	160
Dimensions without condenser (w/d/h)	156 x 236 x 196 mm