

## Unimax 2010

P/N: 036130200



The resilient one with a maximum loading capacity of 10 kg

The Unimax 2010 is the optimal model for increased sample throughput. With a usable area of 390 x 340 mm and the possibility of a multi-story structure, space is used effectively.

The orbit of 20 mm and speeds between 20 and 400 rpm guarantee gentle mixing of all samples. This shaker can be operated either continuously or with an integrated timer (max. 120 min).

## Unimax 2010 - Technical Data

|                                       |                                                                                              |
|---------------------------------------|----------------------------------------------------------------------------------------------|
| Overheat protection                   | self-resetting                                                                               |
| Permissible ambient conditions        | 5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity |
| Weight                                | 16 kg                                                                                        |
| Protection class (EN 60529)           | IP 20                                                                                        |
| Timer                                 | yes (max. 120 min)                                                                           |
| Rotation speed range                  | 20 - 400 rpm                                                                                 |
| Power input                           | 115                                                                                          |
| Dimensions (w/d/h)                    | 426 x 435 x 135 mm                                                                           |
| Max. Load                             | 10 kg                                                                                        |
| Operating Mode                        | timer or continuous                                                                          |
| Rotation speed setting                | digital                                                                                      |
| Orbit / Stroke                        | 20 mm                                                                                        |
| Platform size                         | 390 x 340 mm                                                                                 |
| Accessories included                  | non-skid rubber mat                                                                          |
| Angle                                 | -                                                                                            |
| Motion                                | orbital                                                                                      |
| Ambient temperature range (operation) | -                                                                                            |
| Ambient temperature range (storage)   | -                                                                                            |
| Interfaces                            | -                                                                                            |
| Usable area                           | -                                                                                            |
| Installation altitude                 | -                                                                                            |