


Findenser Mini B14 Cone, B14 Socket

P/N: 015610005



- The high performance Findenser™ Mini is for syntheses under reflux with solvent volumes up to 125 ml.
- It can be secured with standard laboratory clamps or with the Findenser™ attachment kit (sold separately).
- Save water, space and costs with Findenser™ Mini!

Findenser Mini B14 Cone, B14 Socket - Technical Data

Weight	670 kg
Ambient temperature range (operation)	0 °C to 24 °C
Length x Diameter	275 mm x 72 mm
Cone / Joint Size	B14
Max. Solvent Load / Volume	125 ml of solvent in a 250 ml flask
Inner or wetted parts	Borosilicate glass
Finned jacket	Anodised aluminium
Plug & Seal	Acetal & HT Silicone
Thermofluid	Water
Max. solvent bpt	155 °C
Oven drying	Max. Temp. 60 °C
Dishwasher	Max. Temp. 50 °C
Autoclaving	Do Not Autoclave