

BR 14 Collapsible-Blade Impeller

P/N: 036300370



Blade and Half-Moon Impeller

- These impellers are recommended for applications which require average speed
- For mixing tasks with little or average viscosity
- With collapsible blade for narrow neck vessels

BR 14 Collapsible-Blade Impeller - Technical Data

| | |
|---------------------|---------------------|
| Blade size | 90 x 10 mm |
| Material | stainless steel V4A |
| Length | 400 mm |
| Shaft diameter | 8 mm |
| Max. rotation speed | 800 rpm |