


PR 31 Ringed Propeller

P/N: 036300410



Propeller-Type Impeller

- These impellers are recommended for applications which require average or high speed
- For mixing tasks with medium or high viscosity
- Excellent mixing properties for homogenization and suspensions
- These models create an axial flow

PR 31 Ringed Propeller - Technical Data

| | |
|---------------------|----------------------------|
| Blade size | 33 mm |
| Material | stainless steel AISI 316Ti |
| Length | 400 mm |
| Shaft diameter | 8 mm |
| Max. rotation speed | 2,000 rpm |