

Hei-VAC Valve Tec

P/N: 591-00160-00



- Recommended for solvents with low or medium boiling points
- Vacuum can be controlled via valve operated or manual vacuum controllers
- A Vacuum valve is necessary, if you want to combine this vacuum pump with the Hei-VAP Expert Control / Hei-VAP Ultimate Control rotary evaporators

Condenser is not in scope of delivery.

Hei-VAC Valve Tec - Technical Data

Weight	10.5 kg
Diaphragm pump	Two-stage
Suction capacity	1000 m³/h
Ultimate vacuum	0.1 mbar
Combinable with a condenser	yes
Power input	80
Dimensions without condenser (w/d/h)	145 x 338 x 185 mm