

Hei-CHILL Peltier

P/N: 591-01910-00-0



Hei-CHILL Peltier Thermostat 600

Hei-PROCESS
Hei-VAP Expert Control Ultimate Control

-
- 0.05
-
- -10 °C ~ 70 °C
- 2.83
- , RS 232 RS 485

Hei-CHILL Peltier -

	15 kg
	-10 - 70 ° C
Temperature measurement accuracy	0.05 ° C
(/ /)	27 x 43 x 43 cm
Cooling capacity	600 W
Max. flow rate	3 L / min
Max. noise level	63 db
Interfaces	Ethernet, RS 232, RS 485
Control panel	5 “ Touchscreen
	>1000 W