

Reax top

P/N: 541-10000-00



The classic with 5 mm shaking orbit

This model is particularly suitable for simple mixing tasks. In short-term operation mode, the shaking motion is triggered by pressure on the test tube tray. This mode is suitable for mixing up solutions, for example.

Persistent cell pellets can be dissolved gently but powerfully in continuous operation. With a speed of up to 3,200 rpm and a shaking orbit of 5 mm, reliable and uniform distribution of solids, and homogenization of highly viscous media are guaranteed.

Reax top - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	2,8 kg
Protection class (EN 60529)	IP 22
Zamanlayıcı	-
Rotasyon hız derecesi	200 - 3,200 rpm
Güç girişi	51
Dimensions (w/d/h)	134 x 172 x 105 mm
Maksimum yükleme	-
Operasyon modu	automatic or continuous
Rotasyon hız ayarları	analog
Yörünge/ strok	5 mm
Platform ölçüsü	-
Aksesuarlar içerir	-
Açı	-
Hareket	vibrating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-