

Reax control

P/N: 541-11000-00



The precise one with speed scaling

All the outstanding features of the Reax Top's "little brother" can be found in this model: a rotational speed of up to 2,500 rpm and a 5 mm shaking orbit.

The electronic speed control and the scaling of the speed display allow precise settings. Applications that require a determined speed are thus shaken more gently. Shear force-sensitive organisms remain particularly protected by the constant speed in the low range.

Reax control - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	2,8 kg
Protection class (EN 60529)	IP 22
Zamanlayıcı	-
Rotasyon hız derecesi	0 - 2,500 rpm
Güç girişi	51
Dimensions (w/d/h)	134 x 172 x 105 mm
Maksimum yükleme	-
Operasyon modu	automatic or continuous
Rotasyon hız ayarları	analog
Yörünge/ strok	5 mm
Platform ölçüsü	-
Aksesuarlar içerir	-
Açı	-
Hareket	vibrating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-