

Rotamax 120

P/N: 544-41200-00



The compact one with a maximum loading capacity of 2 kg

Space-saving and with a load weight of up to 2 kg, this shaker can perfectly master all applications.

The variably adjustable speed control (20-300 rpm) and the orbit of 20 mm guarantee homogeneous distribution. Your applications can be carried out in continuous operation or by means of a timer (max. 120 minutes).

Rotamax 120 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,5 kg
Protection class (EN 60529)	IP 30
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	20 - 300 rpm
Güç girişi	33
Dimensions (w/d/h)	245 x 310 x 125 mm
Maksimum yükleme	2 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	20 mm
Platform ölçüsü	220 x 220 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	-
Hareket	orbital
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-