

## Unimax 1010

P/N: 543-12310-00



The one with temperature-control and a maximum loading capacity of 5 kg

This shaker can be integrated into the modular incubation system, and it is therefore ideal for all temperature-controlled shaking tasks. In addition, it is possible to operate several shakers of this size in one climate chamber. We recommend climate chambers made by Memmert.

With a speed range of 30-500 rpm and an orbit of 10 mm, samples are optimally kept in motion. Time-dependent applications can be performed without any problems by using the timer (max. 999 min).

## Unimax 1010 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	8 kg
Koruma sınıfı DIN EN 60529	IP 40
Zamanlayıcı	Yes
Rotasyon hız derecesi	30 - 500 rpm
Güç girişi	50 W
Boyutlar (G/ D/ Y)	320 x 375 x 125 mm
Maksimum yükleme	5 kg
Operasyon modu	timer or continuous
Rotasyon hız ayarları	digital
Yörünge/ strok	10 mm
Platform ölçüsü	290 x 258 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	-
Hareket	orbital