

Unimax 2010

P/N: 542-10020-00



The resilient one with a maximum loading capacity of 10 kg

The Unimax 2010 is the optimal model for increased sample throughput. With a usable area of 390 x 340 mm and the possibility of a multi-story structure, space is used effectively.

The orbit of 20 mm and speeds between 20 and 400 rpm guarantee gentle mixing of all samples. This shaker can be operated either continuously or with an integrated timer (max. 120 min).

Unimax 2010 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	16 kg
Protection class (EN 60529)	IP 20
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	20 - 400 rpm
Güç girişi	115
Dimensions (w/d/h)	426 x 435 x 135 mm
Maksimum yükleme	10 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	digital
Yörünge/ strok	20 mm
Platform ölçüsü	390 x 340 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	-
Hareket	orbital
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-