

Duomax 1030 (5° tilt angle)

P/N: 543-32205-00



Containers are safely fixed on the anti-slip rubber mat even at high speeds. The speed is variably adjustable between 2 and 50 rpm. The shaker can be operated in continuous operation or with a timer (max. 120 minutes).

A tilt angle of 5° is recommended for applications for which a gentle movement amplitude is suitable.

Additionally, the Duomax 1030 can be integrated into the modular incubator model.

Duomax 1030 (5° tilt angle) - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	8 kg
Protection class (EN 60529)	IP 20
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	2 - 50 rpm
Güç girişi	115
Dimensions (w/d/h)	320 x 375 x 185 mm
Maksimum yükleme	5 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	-
Platform ölçüsü	290 x 258 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	5°
Hareket	rocking
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-