

Promax 2020

P/N: 542-20020-00



The resilient one with load capacities of 10 kg

With a usable area of 390 x 340 mm, the large Promax is particularly suitable when the use of large or multiple vibrating funnels is required. In conjunction with the appropriate [assembly parts](#), the device can be expanded by a second floor.

This shaker has an integrated timer for up to 120 minutes. An acoustic warning signals the end of the set time and the shaking process is also stopped. A stroke of 20 mm combined with a speed of 20 – 400 rpm guarantees clean phase separation and efficient extraction results.

Promax 2020 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	16 kg
Protection class (EN 60529)	IP 20
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	20 - 400 rpm
Güç girişi	115
Dimensions (w/d/h)	426 x 435 x 125 mm
Maksimum yükleme	10 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	digital
Yörünge/ strok	20 mm
Platform ölçüsü	390 x 340 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	-
Hareket	reciprocating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-