

## Titramax 101

P/N: 544-11300-00



The powerful with a 3 mm shaking orbit

With a capacity of four multi-well plates, the Titramax 101 ensures more intensive mixing results. With a 3 mm orbit, this model is particularly suitable for processes in which solids must be distributed or dissolved in the cavities.

The speed setting of 150-1,350 rpm guarantees a continuously variable adjustment. Use the integrated timer with a runtime of up to 120 minutes for defined processes. The shaker can of course also be used in continuous operation.

## Titramax 101 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,5 kg
Protection class (EN 60529)	IP 30
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	150 - 1,350 rpm
Güç girişi	31
Dimensions (w/d/h)	245 x 310 x 125 mm
Maksimum yükleme	2 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	3 mm
Platform ölçüsü	220 x 220 mm
Aksesuarlar içerir	space for 4 microtier plates
Açı	-
Hareket	vibrating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-