

Titramax 100

P/N: 544-11200-00



The gentle with a 1.5 mm shaking orbit

This shaker is particularly suitable for processes in which cell cultures must be kept in gentle motion.

The nutrients are excellently distributed. This is guaranteed by the 1.5 mm shaker orbit and a speed range of 150-1,350 rpm.

The integrated timer can be set for up to 120 minutes; an acoustic signal sounds when the interval has expired.

Titramax 100 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,5 kg
Protection class (EN 60529)	IP 30
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	150 - 1,350 rpm
Güç girişi	31
Dimensions (w/d/h)	245 x 310 x 125 mm
Maksimum yükleme	2 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	1,5 mm
Platform ölçüsü	220 x 220 mm
Aksesuarlar içerir	space for 4 microtier plates
Açı	-
Hareket	vibrating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-