

Titramax 1000

P/N: 544-12200-00



The one with temperature-control

The largest model of the Titramax series contains all the features of the Titramax 100: a shaking orbit of 1.5 mm, speeds of up to 1,350 rpm, and a timer.

With a storage capacity of up to six multi-well plates, a higher number of samples can be processed. The option to integrate the Titramax 1000 into the modular incubation system also allows temperature-controlled processes to be carried out without problems.

With this model, easy loading and removal of the plates can be carried out comfortably due to the recess and the spacings in the rubber mat.

Titramax 1000 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,5 kg
Koruma sınıfı DIN EN 60529	IP 30
Zamanlayıcı	Yes
Rotasyon hız derecesi	150 - 1,350 rpm
Güç girişi	31 W
Boyutlar (G/ D/ Y)	320 x 375 x 125 mm
Maksimum yükleme	5 kg
Operasyon modu	timer or continuous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	1,5 mm
Platform ölçüsü	290 x 258 mm
Aksesuarlar içerir	space for 6 microtiter plates
Açı	-
Hareket	vibrating