

Vibramax 100

P/N: 544-21200-00



The compact one with a load capacity of 2 kg

Due to the 3 mm orbit and the speed of 150 - 1,350 rpm, this circular vibrating shaker is ideal for suspending, for solubility studies and leaching studies.

The optional [tension roller attachment](#) offers the option of using different container sizes. In combination with the non-slip rubber mat, the vessels are shaken safely even at high speeds. In addition, a timer with up to 120 minutes is integrated, which triggers an acoustic signal after expiration.

Vibramax 100 - Teknik veri

Aşırı ısınma koruması	self resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,5 kg
Protection class (EN 60529)	IP 30
Zamanlayıcı	yes (max. 120 min)
Rotasyon hız derecesi	150 - 1,350 rpm
Güç girişi	31
Dimensions (w/d/h)	245 x 310 x 125 mm
Maksimum yükleme	2 kg
Operasyon modu	timer or continous
Rotasyon hız ayarları	electronic control
Yörünge/ strok	3 mm
Platform ölçüsü	220 x 220 mm
Aksesuarlar içerir	non-skid rubber mat
Açı	-
Hareket	vibrating
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-