

Reax 2

P/N: 541-21001-00



The compact one for small volumes

The Reax 2 with the optional adapter for test tubes is designed for simultaneous mixing of up to 20 vessels with diameters of 10-18 mm. This corresponds to the size of classic test tubes or 15 ml Falcon tubes. Alternatively, you can insert your complete rack into the universal adapter (width: 230 mm, height min./max.: 45 mm/160 mm) of the overhead shaker.

Equipped with a solid metal frame, its design stands for extreme load capacity and reliability. Guaranteed safety during unattended continuous operation.

Reax 2 - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,2 kg
Protection class (EN 60529)	IP 21
Zamanlayıcı	-
Rotasyon hız derecesi	20 - 100 rpm
Güç girişi	27
Dimensions (w/d/h)	510 x 180 x 235 mm
Maksimum yükleme	1 kg
Operasyon modu	-
Rotasyon hız ayarları	analog
Yörünge/ strok	-
Platform ölçüsü	-
Aksesuarlar içerir	-
Açı	-
Hareket	overhead
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-