

## Reax 20 - for 8 bottles (2-32 rpm)

P/N: 541-20008-03



Overhead shaker for 8 Bottles / 2 to 32 rpm (picture similar)

- Set and continuously adjust the variable speed on the analog control knob from 2 to 32 rpm – ideal for waste water analysis
- Easy and fast attachments of bottles
- Other vessels with 160 - 270-mm height and max. 136-mm diameter can be used with attachments

## Reax 20 - for 8 bottles (2-32 rpm) - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	33 kg
Koruma sınıfı DIN EN 60529	IP 21
Zamanlayıcı	-
Rotasyon hız derecesi	2 - 32 rpm
Güç girişi	280 W
Boyutlar (G/ D/ Y)	770 x 520 x 465 mm
Maksimum yükleme	30 kg
Operasyon modu	-
Rotasyon hız ayarları	electronic control
Yörünge/ strok	-
Platform ölçüsü	-
Aksesuarlar içerir	Tension plate for cap dia. up to 77 mm
Açı	-
Hareket	overhead