

## Reax 20 - for 8 bottles (1-16 rpm)

P/N: 541-20008-00



The resilient one with a capacity of up to 8 bottles (picture similar)

The Reax 20 is an ideal shaker for processing and storing samples that tend to sediment. The fields of application range from environmental technology to paint processing and quality assurance in the materials industry.

The maximum permissible weight is 30 kg, so that samples with a high solid content can also be properly mixed. Since this device sets great forces in motion, it is equipped with an emergency stop switch for your safety.

This overhead shaker can be ordered with a capacity of [4](#), [8](#) or [12](#) bottles. Also suitable for mixing cylinders or wide-neck bottles.

## Reax 20 - for 8 bottles (1-16 rpm) - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	28 kg
Protection class (EN 60529)	IP 21
Zamanlayıcı	-
Rotasyon hız derecesi	1 - 16 rpm
Güç girişi	280
Dimensions (w/d/h)	770 x 520 x 465 mm
Maksimum yükleme	30 kg
Operasyon modu	-
Rotasyon hız ayarları	electronic control
Yörünge/ strok	-
Platform ölçüsü	-
Aksesuarlar içerir	Tension plate for cap dia. up to 77 mm
Açı	-
Hareket	overhead
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-