

Kryo 30 liquid (10 l)

P/N: 569-00611-00



Non-flammable thermal fluid for heating and low-temperature range.
10 l thermal fluid for Hei-CHILL Chillers

Kryo 30 liquid (10 l) - Teknik veri

Sıcaklık dercesi	-30 – +90 °C
Chemical characterization	Monoethylene glycol/water mixture containing inhibitors
Color	Light yellow
Viscosity, kinematic at 20 °C	4 mm ² /s
Destiny at 20 °C	1.082 kg/m ³
Pour point	< -70 °C
Flash point	ca. 120 °C
Water solubility	Completely soluble
Material incompatibility	Not known incompatibility with other materials