

## Hei-TORQUE Value 100

P/N: 501-61010-00



Basic model

The reliable solution for all standard applications that do not need documentation

These stirrers are ideal for standard stirring tasks. They are designed to mix and disperse media that require non-reproducible results in high-viscosity applications

- Indication of torque tendencies to monitor viscosity changes
- Digital 2.4" display for ease of operation
- Constant speed even under changing loads
- Safe start and stop of operation via slide touch panel to avoid unintended stirring
- Newest motor generation for maximum power at minimum noise level

## Hei-TORQUE Value 100 - Teknik veri

Güç derecelendirme motor girişi	90 W
Güç derecelendirme motor çıkışı	50 W
Vitesin numarası	1
Rotasyon hız göstergesi	digital monochrom 2.4"
Hız kontrol	electronic
Tork, maksimum	100 Ncm
Tork göstergesi	symbol
Aşırı ısınma koruması	automatic cut-out
Motor koruması	temperature control software
Maks. Viskozite	60.000 mPa s
Maks. karıştırma kapasitesi (H2O)	50 l
Analog/ dijital arayüz	-
Kabul edilebilir durum	continuous operation
Sayaç/ Zamanlayıcı	-
Şaft çapı , maks.	10,5 mm
İzin verilen ortam koşulları	5 – 31 °C at 80 % rel. humidity 32 – 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	4,4 kg
Koruma sınıfı DIN EN 60529	IP 54
Rotasyon hız derecesi	10 - 2,000 rpm
Stay bar ölçüsü (dia. x l)	13 x 160 mm
Boyutlar (G/ D/ Y)	86 x 247 x 350 mm
Rotasyon yön değişimi	-
Boyunca- şaft dizayn	1