

Reax 2 inc. universal adaptor

P/N: 541-21009-00



- Quickly clamps in 2 vessels of any size and accepts load capacities up to 1 kg
- Speed continuously adjustable from 20 to 100 rpm
- Universal adaptor allows for the use of various vessels of 50 - 160-mm height
- In addition, an optional adaptor is available for the use of 20 test tubes at a time

Reax 2 incl. universal adaptor - Teknik veri

Aşırı ısınma koruması	self-resetting
İzin verilen ortam koşulları	5 - 31 °C at 80 % rel. humidity 32 - 40 °C decreasing linearly up to max. 50 % rel. humidity
Ağırlık	5,2 kg
Protection class (EN 60529)	IP 21
Zamanlayıcı	-
Rotasyon hız derecesi	20 - 100 rpm
Güç girişi	27
Dimensions (w/d/h)	510 x 180 x 235 mm
Maksimum yükleme	1 kg
Operasyon modu	-
Rotasyon hız ayarları	analog
Yörünge/ strok	-
Platform ölçüsü	G 230 x Y 198 mm
Aksesuarlar içerir	universal adapter
Açı	-
Hareket	overhead
Ortam sıcaklık dercesi (operasyon)	-
Ortam sıcaklık dercesi (depolanan)	-
Interfaces	-
Usable area	-
Installation altitude	-