

Findenser™ NS19 Cone, NS19 Socket

P/N: 505-81720-00



- The high performance Findenser™ as standard size (400 mm) is for syntheses under reflux with solvent volumes up to 1 liter.
- It can be secured with standard laboratory clamps or with the Findenser™ attachment kit (sold separately).
- Save water, space and costs with Findenser™!

Findenser™ NS19 Cone, NS19 Socket - Teknik veri

Ağırlık	1.200 g
Ortam sıcaklık dercesi (operasyon)	0 °C to 24 °C
Length x Diameter	400 mm x 72 mm
Cone / Joint Size	B19
Max. Solvent Load / Volume	1 litre of solvent in a 2 litre flask
Inner or wetted parts	Borosilicate glass
Finned jacket	Anodised aluminium
Plug & Seal	Acetal & HT Silicone
Thermofluid	Water
Max solvent bpt	155 °C
Oven drying	Max. Temp. 60 °C
Dishwasher	Max. Temp. 50 °C
Autoclaving	Do Not Autoclave