

Findenser™ NS29 Cone, NS24 Socket

P/N: 505-81710-00



- The high performance Findenser™ as standard size (400 mm) is for syntheses under reflux with solvent volumes up to 1 liter.
- It can be secured with standard laboratory clamps or with the Findenser™ attachment kit (sold separately).
- Save water, space and costs with Findenser™!

Findenser™ NS29 Cone, NS24 Socket - Teknik veri

| | |
|------------------------------------|---------------------------------------|
| Ağırlık | 1.200 g |
| Ortam sıcaklık dercesi (operasyon) | 0 °C to 24 °C |
| Length x Diameter | 400 mm x 72 mm |
| Cone / Joint Size | B29 |
| Max. Solvent Load / Volume | 1 litre of solvent in a 2 litre flask |
| Inner or wetted parts | Borosilicate glass |
| Finned jacket | Anodised aluminium |
| Plug & Seal | Acetal & HT Silicone |
| Thermofluid | Water |
| Max solvent bpt | 155 °C |
| Oven drying | Max. Temp. 60 °C |
| Dishwasher | Max. Temp. 50 °C |
| Autoclaving | Do Not Autoclave |