


Findenser™ Mini NS14 Cone, NS14 Socket

P/N: 505-81820-00



- The high performance Findenser™ Mini is for syntheses under reflux with solvent volumes up to 125 ml.
- It can be secured with standard laboratory clamps or with the Findenser™ attachment kit (sold separately).
- Save water, space and costs with Findenser™ Mini!

Findenser™ Mini NS14 Cone, NS14 Socket - Teknik veri

Ağırlık	670 kg
Ortam sıcaklık dercesi (operasyon)	0 °C to 24 °C
Length x Diameter	275 mm x 72 mm
Cone / Joint Size	B14
Max. Solvent Load / Volume	125 ml of solvent in a 250 ml flask
Inner or wetted parts	Borosilicate glass
Finned jacket	Anodised aluminium
Plug & Seal	Acetal & HT Silicone
Thermofluid	Water
Max solvent bpt	155 °C
Oven drying	Max. Temp. 60 °C
Dishwasher	Max. Temp. 50 °C
Autoclaving	Do Not Autoclave