

Hei-VAC Valve Control

P/N: 591-00130-00



- Suction capacity for up to three rotary evaporators at the same time
- Vacuum can be controlled via valve operated or manual vacuum controllers
- Depending on your application you can switch on and off the vacuum pump via switchbox
- A Vacuum valve is necessary, if you want to combine this vacuum pump with the Hei-VAP Expert Control / Hei-VAP Ultimate Control rotary evaporators

Condenser is not in scope of delivery.

Hei-VAC Valve Control - Teknik veri

Ağırlık	12,8 kg
Diyagram pompa	2-aşama
Emme kapasitesi	2 m ³ /h
En üst düzey vakum	7 mbar
Kondenser ile birleştirilebilir	Evet
Güç girişi	180
Dimensions without condenser (w/d/h)	195 x 245 x 310 mm