

Hei-VAC Valve Tec

P/N: 591-00160-00



- Recommended for solvents with low or medium boiling points
- Vacuum can be controlled via valve operated or manual vacuum controllers
- A Vacuum valve is necessary, if you want to combine this vacuum pump with the Hei-VAP Expert Control / Hei-VAP Ultimate Control rotary evaporators

Condenser is not in scope of delivery.

Hei-VAC Valve Tec - Teknik veri

Ağırlık	6,5 kg
Diyagram pompa	2-aşama
Emme kapasitesi	0,75 m ³ /h
En üst düzey vakum	12 mbar
Kondenser ile birleştirilebilir	Evet
Güç girişi	80
Dimensions without condenser (w/d/h)	145 x 338 x 185 mm