

Hei-VAC Vario Station

P/N: 591-00142-00



Fully controllable stand-alone pumping unit including vacuum controller.
For Hei-VAP Core models

- Precise vacuum control avoids bumping and eliminates foaming of your evaporation solution
- Due to automatic vacuum supply to the process parameters you have more time for important laboratory activities

Condenser is not in scope of delivery.

Hei-VAC Vario Station - Teknik veri

Ağırlık	6 kg
Diyagram pompa	3-aşama
Emme kapasitesi	1,7 m ³ /h
En üst düzey vakum	5 mbar
Kondenser ile birleştirilebilir	Evet
Güç girişi	160
Dimensions without condenser (w/d/h)	294 x 192 x 299 mm