

Hei-VAC Vario Control

P/N: 591-00141-00



For Hei-VAP Expert Control / Hei-VAP Ultimate Control

Even if the gas ballast valve is open an excellent ultimate vacuum is reached when working with easily condensable vapors. This makes distilling high boiling point solvents such as DMF or DMSO possible at low bath temperatures.

Condenser is not in scope of delivery.

Rotavac Vario Control - Technical Data

Ağırlık	5,4 kg
Diyagram pompa	3-aşama
Emme kapasitesi	1,7 m ³ /h
En üst düzey vakum	2 mbar
Kondenser ile birleştirilebilir	Evet
Güç girişi	160
Dimensions without condenser (w/d/h)	167 x 236 x 196 mm