


BR 10 Cross-Blade Impeller

P/N: 509-10000-00



Due to its blade dimensions, the cross-blade impeller BR 10 is ideally suited for solutions, emulsions and suspensions in vessels with the following characteristics:

- ideal vessel diameter: 100 mm
- ideal ground clearance: 5 mm
- ideal filling level: 100 mm
- ideal volume: 800 ml

The cross-blade impeller BR 10 has a total length of 400 mm and is made of stainless steel. Its stirring shaft has a diameter of 8 mm.

BR 10 Cross-Blade Impeller - Teknik veri

Blade ölçüsü	50 x 12 mm
Materyal	stainless steel V4A
Uzunluk	400 mm
Şaft çapı	8 mm
Mak. Rotasyon hızı	2.000 rpm