

BR 14 Collapsible-Blade Impeller

P/N: 509-14000-00



Especially when using narrow-necked vessels, e.g. round-bottom flasks, an stirring tool with collapsible blades is used. Just like our collapsible-blade impeller BR 14.

This blade impeller is ideal for dissolving and homogenizing in vessels with the following characteristics:

- ideal vessel diameter: 180 mm
- ideal ground clearance: 9 mm
- ideal filling level: 180 mm
- ideal volume: 4,600 ml

The collapsible-blade impeller BR 14 has a total length of 400 mm and is made of stainless steel. Its stirring shaft has a diameter of 8 mm.

BR 14 Collapsible-Blade Impeller - Teknik veri

Blade ölçüsü	90 x 10 mm
Materyal	stainless steel V4A
Uzunluk	400 mm
Şaft çapı	8 mm
Mak. Rotasyon hızı	800 rpm