

PR 30 Pitched-Blade Impeller

P/N: 509-30000-00



Due to its blade dimensions, the pitched-blade impeller PR 30 is ideal for solutions, emulsions and suspensions in vessels with the following characteristics:

- ideal vessel diameter: 145- 290 mm
- ideal ground clearance: 29- 58 mm
- ideal filling level: 145- 290 mm
- ideal volume: 2,400- 19,200 ml

The pitched-blade impeller PR 30 has a total length of 400 mm and is made of stainless steel. Its stirring shaft has a diameter of 8 mm.

PR 30 Pitched-Blade Impeller - Teknik veri

| | |
|--------------------|---------------------|
| Blade ölçüsü | 58 mm |
| Materyal | stainless steel V4A |
| Uzunluk | 400 mm |
| Şaft çapı | 8 mm |
| Mak. Rotasyon hızı | 2.000 rpm |