


PR 31 Ringed Pitched-Blade Impeller (33 mm)

P/N: 509-31000-00



Due to its blade dimensions, the ringed pitched-blade impeller PR 31 is ideal for solutions, emulsions and suspensions in vessels with the following characteristics:

- ideal vessel diameter: 83- 165 mm
- ideal ground clearance: 17- 33 mm
- ideal filling level: 83- 165 mm
- ideal volume: 400- 3,500 ml

The ringed pitched-blade impeller PR 31 has a total length of 400 mm and is made of stainless steel. Its stirring shaft has a diameter of 8 mm.

PR 31 Ringed Pitched-Blade Impeller (33 mm) - Teknik veri

Blade ölçüsü	33 mm
Materyal	stainless steel V4A
Uzunluk	400 mm
Şaft çapı	8 mm
Mak. Rotasyon hızı	2.000 rpm