

## PR 32 Ringed Propeller

P/N: 509-32000-00



### Propeller-Type Impeller

- These impellers are recommended for applications which require average or high speed
- For mixing tasks with medium or high viscosity
- Excellent mixing properties for homogenization and suspensions
- These models create an axial flow

## PR 32 Ringed Propeller - Teknik veri

Blade ölçüsü	45 mm
Materyal	stainless steel V4A
Uzunluk	400 mm
Şaft çapı	8 mm
Mak. Rotasyon hızı	2.000 rpm