

## TR 20 Radial-Flow Impeller (28 mm)

P/N: 509-20000-00



Due to the jagged shape of the stirrer blade, the radial flow impeller TR 20 is ideally suited for the construction of emulsions and the gassing or foaming of liquids in vessels with the following characteristics:

- ideal vessel diameter: 80- 140 mm
- ideal ground clearance: 14- 28 mm
- ideal filling level: 80- 140 mm
- ideal volume: 400- 2,200 ml

The radial-flow impeller TR 20 has a total length of 400 mm and is made of stainless steel. Its stirrer shaft has a diameter of 8 mm.

## TR 20 Radial-Flow Impeller (28 mm) - Teknik veri

Blade ölçüsü	28 mm
Materyal	stainless steel V4A
Uzunluk	400 mm
Şaft çapı	8 mm
Mak. Rotasyon hızı	2.000 rpm