

VISCO JET® - 60 mm
P/N: 509-16060-00



- Reduce your process times significantly while performing the best mixing results ever
- One system for literally all mixing tasks for low to high-viscosity media
- The turbulent flow which is created by a special cone principle even at low speeds is unique to the VISCO JET®
- Even with high-viscosity media and gels which naturally do not mix when common impellers are used you will observe an immediate flow through the entire beaker
- This technology allows for de-gassing of gels while preventing air intake and foaming
- A stirrer shaft is included as standard.

VISCO JET® - 60 mm - Teknik veri

Blade ölçüsü	60 mm
Materyal	stainless steel V4A
Uzunluk	500 mm
Şaft çapı	10 mm
Vesel çapı için	80 - 150 mm
Rotasyon hız derecesi	200 - 800 rpm